

Ambu® Aura™ family

- Quality that runs in the family



Ambu® AuraGain™

The AuraGain is Ambu's new laryngeal mask, satisfying 3 fundamental airway management needs by integrating gastric access and intubation capability in an anatomically curved single-use device that facilitates rapid establishment of a safe airway.

The original anatomical curve is pre-formed to follow the anatomy of the human airway, and the soft rounded curve of the AuraGain ensures rapid placement and guarantees long-term performance.

The thin and soft cuff of the AuraGain is designed to deliver high seal pressures - documented up to 40 cm H₂O.*

The integrated gastric access channel is designed with a low friction inner surface to facilitate easy placement of a gastric tube. Introduce a gastric tube through the device and into the stomach of the patient to enable active and passive management of gastric content, and prevent gastric insufflation.

The AuraGain provides the added safety feature of intubation capability. This means that in case of an unexpected difficult airway, or a "Cannot Intubate – Cannot mask Ventilate" (CI-CV) situation, where the end-game is to intubate the patient, AuraGain can be used as conduit for direct endotracheal intubation assisted by a flexible scope (such as the Ambu® aScope™ 4).

*data on file

Technical specifications and ordering information

| Item no. | Size | Patient weight | Max. inflation volume cuff | Max. intra-cuff pressure | Max. gastric tube | Max. ETT tube | Connector | Material |
|-----------|------|----------------|----------------------------|--------------------------|-------------------|---------------|-----------|--------------------|
| 408100000 | 1 | <5 kg | 4 ml | 60 cm H ₂ O | 6 FR | 3.5 | 15 mm | Phthalate-free PVC |
| 408150000 | 1½ | 5-10 kg | 7 ml | 60 cm H ₂ O | 8 FR | 4.0 | 15 mm | Phthalate-free PVC |
| 408200000 | 2 | 10-20 kg | 10 ml | 60 cm H ₂ O | 10 FR | 5.0 | 15 mm | Phthalate-free PVC |
| 408250000 | 2½ | 20-30 kg | 14 ml | 60 cm H ₂ O | 10 FR | 5.5 | 15 mm | Phthalate-free PVC |
| 408300000 | 3 | 30-50 kg | 20 ml | 60 cm H ₂ O | 16 FR | 6.5 | 15 mm | Phthalate-free PVC |
| 408400000 | 4 | 50-70 kg | 30 ml | 60 cm H ₂ O | 16 FR | 7.5 | 15 mm | Phthalate-free PVC |
| 408500000 | 5 | 70-100 kg | 40 ml | 60 cm H ₂ O | 16 FR | 8.0 | 15 mm | Phthalate-free PVC |
| 408600000 | 6 | >100 kg | 50 ml | 60 cm H ₂ O | 16 FR | 8.0 | 15 mm | Phthalate-free PVC |



Ambu® Aura-i™

Ambu® Aura-i™ is a disposable laryngeal mask designed for everyday use. Pre-formed to follow the anatomy of the human airway, the soft rounded curve of the Aura-i ensures fast and easy placement and guarantees long-term performance.

Aura-i is the preferred everyday routine device for establishing an airway where an ET-tube is deemed unnecessary. At the same time it holds a place on the difficult airway cart.

In case of a "Cannot Intubate – Cannot mask Ventilate" (CI-CV) situation, where the end-game is to intubate the patient, Aura-i can be used as conduit for direct endotracheal intubation assisted by a flexible scope (such as the Ambu® aScope™ 3).

Technical specifications and ordering information

| Item no. | Size | Patient weight | Max. inflation volume cuff | Max. intra-cuff pressure | Max. ETT tube | Connector | Material |
|-------------|------|----------------|----------------------------|--------------------------|---------------|-----------|--------------------|
| 329 100 000 | 1 | <5 kg | 4 ml | 60 cm H ₂ O | 3.5 | 15 mm | Phthalate-free PVC |
| 329 150 000 | 1½ | 5-10 kg | 7 ml | 60 cm H ₂ O | 4.0 | 15 mm | Phthalate-free PVC |
| 329 200 000 | 2 | 10-20 kg | 10 ml | 60 cm H ₂ O | 5.0 | 15 mm | Phthalate-free PVC |
| 329 250 000 | 2½ | 20-30 kg | 14 ml | 60 cm H ₂ O | 5.5 | 15 mm | Phthalate-free PVC |
| 329 300 000 | 3 | 30-50 kg | 20 ml | 60 cm H ₂ O | 6.5 | 15 mm | Phthalate-free PVC |
| 329 400 000 | 4 | 50-70 kg | 30 ml | 60 cm H ₂ O | 7.5 | 15 mm | Phthalate-free PVC |
| 329 500 000 | 5 | 70-100 kg | 40 ml | 60 cm H ₂ O | 8.0 | 15 mm | Phthalate-free PVC |
| 329 600 000 | 6 | >100 kg | 50 ml | 60 cm H ₂ O | 8.0 | 15 mm | Phthalate-free PVC |



Ambu® AuraOnce™

Ambu AuraOnce is a disposable laryngeal mask for everyday use. AuraOnce has a special curve that carefully replicates natural human anatomy. Molded directly into the tube, the unique shape makes correct insertion fast and easy.

AuraOnce's curve also ensures that the patient's head remains in a natural supine position when the mask is in use without extra stress on the upper jaw. Moreover, AuraOnce is specially designed to give the airway tube the flexibility needed to adapt to individual anatomical variances and a wide range of head positions.

AuraOnce has a D-shaped airway tube to give a firm and ergonomical grip during insertion. In short, positioning is fast and accurate, with low cuff pressure.

Technical specifications and ordering information

| Item no. | Size | Patient weight | Max. inflation volume cuff | Max. intra-cuff pressure | Connector | Material |
|-------------|------|----------------|----------------------------|--------------------------|-----------|--------------------|
| 321 100 000 | 1 | <5 kg | 4 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 321 150 000 | 1½ | 5-10 kg | 7 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 321 200 000 | 2 | 10-20 kg | 10 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 321 250 000 | 2½ | 20-30 kg | 14 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 321 300 000 | 3 | 30-50 kg | 20 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 321 400 000 | 4 | 50-70 kg | 30 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 321 500 000 | 5 | 70-100 kg | 40 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 321 600 000 | 6 | >100 kg | 50 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |



Ambu® AuraStraight™

Ambu AuraStraight is a disposable laryngeal mask designed for use both in anesthesia and emergency medicine. It provides users with an even greater choice and is an alternative to the already popular AuraOnce curved laryngeal mask.

Unlike traditional laryngeal masks, AuraStraight features a unique one-piece-mold design, which ensures not only consistency, but also that there are no rough edges on the airway tube or cuff.

Furthermore, the AuraStraight provides an extra soft, thin cuff. Thus, the seal more readily conforms to the shape of the airway with the minimal amount of pressure.

Technical specifications and ordering information

| Item no. | Size | Patient weight | Max. inflation volume cuff | Max. intra-cuff pressure | Connector | Material |
|-------------|------|----------------|----------------------------|--------------------------|-----------|--------------------|
| 324 100 000 | 1 | <5 kg | 4 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 324 150 000 | 1½ | 5-10 kg | 7 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 324 200 000 | 2 | 10-20 kg | 10 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 324 250 000 | 2½ | 20-30 kg | 14 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 324 300 000 | 3 | 30-50 kg | 20 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 324 400 000 | 4 | 50-70 kg | 30 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 324 500 000 | 5 | 70-100 kg | 40 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 324 600 000 | 6 | >100 kg | 50 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |



Ambu® AuraFlex™

Ambu AuraFlex is a disposable, flexible laryngeal mask which is specially designed for ENT, ophthalmic, dental and other head and neck surgeries.

AuraFlex has a long, wire-reinforced, and flexible airway tube that allows positioning away from the surgical field without loss of seal. The kink-free wire-reinforced airway tube eliminates the risk of airway tube occlusion.

Furthermore, the integrated pilot tube makes handling and insertion easy and reduces the risk of cutting the pilot tube during surgery.

Technical specifications and ordering information

| Item no. | Size | Patient weight | Max. inflation volume cuff | Max. intra-cuff pressure | Connector | Material |
|-------------|------|----------------|----------------------------|--------------------------|-----------|--------------------|
| 327 200 000 | 2 | 10-20 kg | 10 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 327 250 000 | 2½ | 20-30 kg | 14 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 327 300 000 | 3 | 30-50 kg | 20 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 327 400 000 | 4 | 50-70 kg | 30 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 327 500 000 | 5 | 70-100 kg | 40 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |
| 327 600 000 | 6 | >100 kg | 50 ml | 60 cm H ₂ O | 15 mm | Phthalate-free PVC |



Ambu® Aura40™

Ambu Aura40 is a reusable laryngeal mask for everyday use. Aura40 has a special curve that carefully replicates natural human anatomy. Molded directly into the tube, the unique shape makes correct insertion fast and easy.

Aura40's curve also ensures that the patient's head remains in a more natural position when the mask is in use without extra stress on the upper jaw. Moreover, Aura40 is specially designed to give the airway tube the flexibility needed to adapt to individual anatomical variances and a wide range of head positions.

Technical specifications and ordering information

| Item no. | Size | Patient weight | Max. inflation volume cuff | Max. intra-cuff pressure | Connector | Material |
|-------------|------|----------------|----------------------------|--------------------------|-----------|----------|
| 326 100 000 | 1 | <5 kg | 4 ml | 60 cm H ₂ O | 15 mm | Silicone |
| 326 150 000 | 1½ | 5-10 kg | 7 ml | 60 cm H ₂ O | 15 mm | Silicone |
| 326 200 000 | 2 | 10-20 kg | 10 ml | 60 cm H ₂ O | 15 mm | Silicone |
| 326 250 000 | 2½ | 20-30 kg | 14 ml | 60 cm H ₂ O | 15 mm | Silicone |
| 326 300 000 | 3 | 30-50 kg | 20 ml | 60 cm H ₂ O | 15 mm | Silicone |
| 326 400 000 | 4 | 50-70 kg | 30 ml | 60 cm H ₂ O | 15 mm | Silicone |
| 326 500 000 | 5 | 70-100 kg | 40 ml | 60 cm H ₂ O | 15 mm | Silicone |
| 326 600 000 | 6 | >100 kg | 50 ml | 60 cm H ₂ O | 15 mm | Silicone |



Ambu® Aura™ family

- tried, tested, trusted

The Aura Family is a complete range of easy-to-use SGAs offering high efficiency and increased patient safety for all airway management needs.

The Ambu® legacy

Since 1937, breakthrough ideas have fuelled Ambu's work on bringing efficient healthcare solutions to life. The manifestations of our efforts range from early inventions like the Ambu bag and the legendary BlueSensor™ electrodes to the world's first single-use flexible video scope – the aScope™.

The Ambu® Aura™

For more than a decade, Ambu has specialised in developing SGAs. Moulded in a single piece, and featuring a novel anatomical curve, AuraOnce™ has delivered easy and fast insertion in millions of cases. Aura-i™ is designed for everyday use, but features a built-in backup plan if needed; intubation capability with your standard ET tubes.

Our most recent innovation - AuraGain™ - integrates gastric access and intubation capability in an anatomically curved single-use laryngeal mask.

Over time, the Aura™ family has expanded to cover a comprehensive range of high quality SGAs for your needs. All members of the Aura family are made in soft materials that perform in routine as well as advanced cases.

Every single one of our masks is designed with patient safety and user friendliness in mind, and our solutions are known for adding new value to our customers. We call it the Ambu Aura.

Innovative design that improves patient care

The Aura Family offers a choice for all needs, ranging from a simple straight device to the most advanced 2nd generation SGA.

Tried, tested and trusted by clinicians

The Aura Family offers a full range of clinically proven SGAs with millions of uses all over the world.

Complete range enabling the right choice

The Aura Family is the most versatile range of SGAs for routine as well as more advanced cases.

ambu.com/AuraFamily



US: Rx only

Ambu A/S • Baltorpbakken 13 • DK - 2750 Ballerup • Denmark

Tel. +45 7225 2000 • Fax +45 7225 2053