# INSTRUCTIONS FOR USE

**Ambu**

**Ambu® Open Cuff Silicone Face Mask**

Reusable Mask

## Contents

## Page

English(Instructions for use) .....	1-4
Česky(Návod k použití) .....	5-8
Dansk(Brugsanvisning) .....	9-12
Deutsch(Bedienungsanleitung) .....	13-16
Ελληνικά(Δηγίες χρήσεως) .....	17-20
Español(Manual de instrucciones) .....	21-24
Eestlane(Kasutusjuhendit) .....	25-28
Suomi(Käyttöohje) .....	29-32
Français(Mode d'emploi) .....	33-36
Magyar(Használati útmutató) .....	37-40
Italiano(Manuale d'uso) .....	41-44
日本語(使用法) .....	45-48
Nederlands(Gebruiksaanwijzing) .....	49-52
Norsk(Brukerveiledning) .....	53-56
Polski(Instrukcja obsługi) .....	57-60
Português(Manual de instruções) .....	61-64
Slovenčina(Návod na použitie) .....	65-68
Svenska(Instruktionshandbok) .....	69-72
Türkçe(Kullanım talimatları) .....	73-76

## Ambu® Open Cuff Silicone Face Mask

### Intended Use:

Open Cuff Silicone Face Mask Reusable Face Mask is intended to serve as an interface for delivering respiratory gases from a resuscitator or anesthesia machine to a patient's airway.

### Indications:

- Any situation in which the patient requires ventilation.
- Patients undergoing surgical procedures.

### Contraindications:

- Facial, esophageal, or gastric surgery.
- Trauma or burns involving the face.

### Precautions:

- Inspect the mask for damage before use.
- Be aware the possible air leakage from the mask happened by improper setup.
- Make sure the mask is properly assembled and securely attached to the equipment before use.
- Monitor the equipment system during use to ensure proper function and detect any issues promptly.
- Store the mask in a clean, dry location between uses, and avoid exposing it to extreme temperatures or other damaging conditions.

### Warnings:

- Intended for use by qualified medical personnel.
- Always test this product for proper function before use on a patient.
- DO NOT USE any hydrocarbon-based substance on any part of this device.
- DO NOT clean the mask with alcohol, chlorine or bleach.
- DO NOT store the mask in direct sunlight.
- DO NOT use a mask that is past its expiration date. This can compromise the integrity and effectiveness of the circuit and lead to patient harm.
- DO NOT use mask that have been stored improperly or exposed to extreme temperatures or other damaging conditions.
- Please report any serious incident that has occurred in relation to the device to the manufacturer and competent authority of the member state.

### Pre-use Check:

- Check the mask and circuit for any visible damage or wear and tear, such as cracks, leaks, or holes.
- Check that all connections and fittings are securely fastened and there are no loose components.
- Ensure that the mask is clean and free of any debris or obstructions.

**Directions for Use:**

- Step 1. Remove plastic packaging from Open Cuff Silicone Face Mask Reusable Face Mask.
- Step 2. Select an appropriately sized mask for the patient's face to ensure a proper seal.
- Step 3. Attach the mask to the resuscitator by securely connecting the mask to the patient valve.
- Step 4. Make sure there are no air leaks in the connections.
- Step 5. If you are using supplemental oxygen, adjust the flow rate as needed to achieve the desired oxygen concentration.
- Step 6. Hold the mask firmly against the patient's face, covering the nose and mouth.
- Step 7. With the mask sealed on the patient's face, squeeze the bag of the resuscitator to deliver a breath.
- Step 8. Continuously monitor the patient's chest rise and fall to assess the effectiveness of ventilation.

**Cleaning and Disinfection:****Preparation for Decontamination:**

Remove all connector adapters from the device. The devices must be reprocessed in a disassembled state, as far as possible.

**Pre-Cleaning:**

Do a manual pre-cleaning, until the devices are visually clean. Submerge the devices in a cleaning solution. Clean the surfaces with a soft bristol brush.

**Cleaning:**

Regarding cleaning/disinfection, rinsing and drying, it is to distinguish between manual and automated reprocessing methods. Preference is to be given to automated reprocessing methods, especially due to the better reproducibility and standardisation, and in personnel protection.

**Automated Cleaning:**

1. Use a washer-disinfector meeting the requirements of the ISO 15883 series.
2. Put the instrument into the machine on a tray. Connect the parts with lumen with the WD by using suitable adapter.
3. The devices in the washer-disinfector are arranged in such a way that there is no rinsing shadow and the water drains off quickly, especially for the lumen devices.
4. Start the program:
  - 4 min pre-washing with cold water (<40°C);
  - emptying
  - 5 min washing with a mild alkaline cleaner at 55°C
  - emptying
  - 3 min neutralising with warm water (>40°C);
  - emptying

- 5 min intermediate rinsing with warm water (>40°C)
- Emptying

*The automated cleaning processes have been validated by using 0.5% neodisher MediClean forte(Dr. Weigert). Note Acc. to ISO 17664 no manual reprocessing methods are required for these devices. If a manual reprocessing method has to be used, please validate it prior to use.*

**Disinfection:**

1. Automated thermal disinfection in washer/disinfector under consideration of national requirements in regards to A0 value (see ISO 15883).
2. A disinfection cycle of 5 min disinfection at 90°C has been validated for the device to achieve an A0 value of > 3000. Here we suggest a disinfection cycle of 5 min disinfection time at 93 °C.

**Drying:****Automated Drying:**

Drying the instrument through drying cycle of washer/disinfector. If needed, additional manual drying can be performed through lint free towel. Insufflate cavities of products by using sterile compressed air.

**Functional Testing, Maintenance:**

1. Visual inspection for cleanliness of the products and reassembling, if required.
2. All instruments should be checked again for dryness.
3. After cleaning and disinfection, a thorough inspection and maintenance ensures that the products are fit for use.
4. Check that the product has no dents, cracks, deformations, scratches, etc.
5. Check all markings on the product for clear visibility.
6. Discard and replace any components as necessary.
7. Do not use the device with following defects: material deformation, cracks on the product, brittle or other change in the material, etc.

**Packaging:**

Pack the products in an appropriate packaging material for sterilization. The packaging material and system refer to ISO 11607.

**Sterilization:**

1. Sterilization of products by applying a fractionated pre-vacuum steam sterilization process (according to EN 285/ISO 17665) under consideration of the respective country requirements.
2. Following sterilization parameters are commonly used: 134 °C, 5 min (for product only contain Silicone, PSF) or 121°C, 30min (other reusable material product).

### 3. Drying time:

For steam sterilization, we recommend a drying time of 20 to 40 minutes. Choose a suitable drying time, depending on the autoclave and load. Refer to the autoclave's instructions for use.

### 4. After sterilization:

- Remove the product from the autoclave.
- Let the product cool down at room temperature for at least 30 minutes. Do not use additional cooling.
- Check that the sterilization wraps or pouches are not damaged.

**Duration of Use:** ≤7 Days, no more than 20 times of reuse.

### Storage:

Storage of sterilized instruments in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.

### Disposal:

Dispose of all materials in accordance with hospital protocols and procedures and all applicable local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

**Specifications:** Silicone

### Symbols :

Caution

Catalogue number

Date of manufacture

Use-by date

Batch code

Medical device

Manufacturer

Distributor

Consult instructions for use or consult electronic instructions for use

Authorized representative in the European Community

Importer

**Rx only** Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician

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Ambu and other trademarks are trademarks of Ambu A/S
GaleMed (Xiamen) Co., Ltd.
Ambu A/S
**Emergo Europe B.V.**