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

Attention! The operating instructions must be read and applied prior to operation.

Isolating transformers for the “pro-cart” mobile equipment carriers

Manufacturer:

 **NORATEL**
Germany AG
Elsenthal 53
DE 94481 Grafenau

Table of content

1.0	Description of the isolating transformer	4
1.1	Important safety instructions for use.....	6
1.2	Warnings.....	7
1.3	Explanation of graphic symbols.....	8
1.4	Intended use	9
1.5	Instructions for use	10
1.6	Operating manual	10
1.7	Product features and designs.....	11
1.8	Power supply cords, jumper cords	11
1.9	Changing the fuses.....	12
1.10	Operating conditions.....	12
1.11	Transport and storage conditions	12
1.12	Cleaning, disinfection.....	13
2.0	Description of the integrated insulation monitor (earth-leakage guard) ELG	14
2.1	Intended use	14
2.2	Product characteristics	14
2.3	Operating manual earth-leakage guard	15
2.4	Technical data of the ELG insulation monitor.....	16
3.0	Elimination of faults.....	17
3.1	Information about electromagnetic compatibility (EMC).....	18
3.2	Guarantee conditions.....	19

1.0 Description of the isolating transformer


This isolating transformer is an optional accessory for all models of the ITD “pro-cart” mobile equipment carriers. It is intended for installation in the chassis and is available in the following performance classes:

Item No.	Drawing No.	Description	Output power [VA]	Input voltage [V] ^{*)}	Output voltage [V] ^{*)}	Power consumption [VA]	Weight [kg]
9-059-130128	ZV.9386.999	pro-cart 660 VA	660	115/230	115/230	690	9,5
9-059-130209	ZV.9575.999	pro-cart 660 VA ELG	660	115/230	115/230	690	9,8
9-059-130129	ZV.9387.999	pro-cart 1200 VA	1200	115/230	115/230	1240	14,2
9-059-130210	ZV.9576.999	pro-cart 1200 VA ELG	1200	115/230	115/230	1240	14,5
9-059-130212	ZV.9581.999	pro-cart 1600 VA	1600	115/230	115/230	1650	17,4
9-059-130213	ZV.9582.999	pro-cart 1600 VA ELG	1600	115/230	115/230	1650	17,7

^{*)}Operating frequency 50/60 Hz

The following illustrations show the essential operating and connection elements of the 660 VA and 1200 VA models:

9-059-130128	ZV.9386.999	pro-cart 660 VA
9-059-130209	ZV.9575.999	pro-cart 660 VA ELG
9-059-130129	ZV.9387.999	pro-cart 1200 VA
9-059-130210	ZV.9576.999	pro-cart 1200 VA ELG


- 1 Appliance outlet 6-fold IEC-F
- 2 RJ-45 connector for connecting the external insulation monitor control and display element
- 3 Interface for external on/off switch
- 4 POAG-bolt pursuant to DIN 42801
- 5 System grounding cable
- 6 Appliance inlet IEC-C14
- 7 Fuse holder for the input side fuse inserts
- 8 Power switch green, illuminated
- 9 Primary voltage selector
- 10 Secondary voltage selector

The following illustrations show the essential operating and connection elements of 1600 VA models:

9-059-130212
9-059-130213

ZV.9581.999
ZV.9582.999

pro-cart 1600 VA
pro-cart 1600 VA ELG

- ① Appliance outlet 6-fold IEC-F
- ② RJ-45 connector for connecting the external insulation monitor control and display element
- ③ Interface for external on/off switch
- ④ POAG-bolt pursuant to DIN 42801
- ⑤ System grounding cable
- ⑥ Appliance inlet IEC-C14
- ⑦ Fuse holder for the input side fuse inserts
- ⑧ Primary voltage selector
- ⑨ Secondary voltage selector
- ⑩ Power switch green, illuminated

1.1 Important safety instructions for use

To avoid harmful effects on the user or patients, for safe operation and to avoid damage to the isolating transformer, please observe the following safety instructions:

- Protect the isolating transformer from direct exposure to moisture and wetness.
- Operate the isolating transformer only in dry areas.
- When transferring to a heated room, ensure that the isolating transformer is not put into operation until the isolating transformer has adjusted to the room temperature (about 30 minutes).
- Keep the isolating transformer away from sources of heat (heating units, solar radiation, fan heaters, etc.).
- Due to the thermal radiation emanating from the isolating transformer at full load, temperature-sensitive devices which are located in the immediate vicinity can be impaired. In this case, please increase the distance between the isolating transformer and the device.
- Make sure that the isolating transformer is only connected to a supply mains with a functional protective conductor connection which complies with the provisions of DIN VDE 0100-710:2012-10 respectively IEC 60364-7-710 2012-10 "Electrical installations in buildings – Part 7-710 Requirements for special installations or rooms – rooms used for medical purposes". If in doubt, contact a specialist electric workshop or an authorised hospital technician.
- Ensure that the connection cables do not show any pinch points, bare areas or other damages. If you notice any damage, the affected cable has to be replaced immediately. For this purpose, please contact the supplier or manufacturer of the equipment cart.
- There are no user-serviceable parts in the isolating transformer. Therefore, the unit may only be opened by authorised personnel.
- The isolating transformer is not intended for use in sterile rooms. No sterilisation may be performed. For cleaning, use only a slightly damp cloth with a small amount of mild detergent so that no moisture can penetrate into the isolating transformer.
- It is possible to disinfect the housing surface. However, it is also important to ensure that only slightly damp cloths are used, and that no moisture or wetness can penetrate into the isolating transformer.
- The isolating transformer fulfils all the requirements for outgoing EMC disturbances. If there is any interference with the connected devices, check the wiring for proper operation. Avoid that signal-carrying lines, for example, are laid in the immediate vicinity, parallel to the power supply order or jumper cords or bundled together. If interference occurs, increase the distance between the cables and/ or the device concerned and the isolating transformer. EM disturbances have no influence on the essential performance of the device.
Maintain a distance of approx. 1.5 meters from the magnet system of MRI systems.
Do not use the equipment trolley base frame as a base for DECT or cell phones.
- The connection of the control and display panel of the insulation monitor is only to be performed on the RJ-45 interface which has been provided on the isolating transformer for this purpose. Do not try to connect other components to this interface.
- **The system has to be switched off and disconnected from the supply mains during all work on the isolating transformer or insulation monitor. Secure the system against being switched on unintentionally.**

1.2 Warnings

RISK OF EXPLOSION: There is a risk of explosion if the isolating transformer is operated or used in the vicinity of highly inflammable anaesthetic gases or similar combustible gases.

ATTENTION

This is a device of Protection Class I.
The safe grounding of the isolating transformer has to be checked repeatedly at regular intervals!

ATTENTION

The isolating transformer must be subjected to a safety inspection (ST) at regular intervals. A cycle of at least 2 years is recommended.

ATTENTION

Never open the housing! When the housing is open, there is a risk of death from electric shock. There are no user-replaceable electronic components inside the isolating transformer. Leave repairs or maintenance exclusively to your authorised dealer or to the manufacturer of the mobile equipment carrier. The isolating transformer may not be exposed to water, water vapour or high humidity.

WARNING

Under no circumstances may the isolating transformer be operated with input voltages other than those specified in section 1.0.

ATTENTION

The isolating transformer requires adequate air convection for cooling purposes. Make sure that the air slots are not obstructed.

ATTENTION

Make sure that only the fuse values stated in section 1.9 are used according to the input voltages. Otherwise, the isolating transformer and the associated consumers may fail.

ATTENTION

Anyone connecting additional devices is a system configurator and is therefore responsible for complying with IEC/EN 60601-1, section 16 for the ME-system.

ATTENTION

Make sure that no system-independent devices are connected which do not correspond to the purpose of the medical electrical system. Verify that the connected devices are medical devices tested in accordance with the applicable standards . A change – such as by replacing, eliminating or supplementing individual devices – leads to the medical electrical system having to be reassessed by the system configurator.

1.3 Explanation of graphic symbols

Symbols located on the isolating transformer:

Manufacturer of the isolating transformer

"1" The isolating transformer is switched on, the rocker is illuminated green
"0" The isolating transformer is switched off, the rocker is not illuminated

The triangle with the exclamation point is intended to alert the user of important functions and mainly to refer to maintenance regulations which can be read in the operating instructions.

Please follow the operating instructions!

Connection for the potential equalisation (POAG)
The additional potential equalisation has the task of reducing the potentials of different metal parts that can be touched at the same time or reducing potential differences that may arise during application between the body, electromedical devices and external conductive parts.



This is not an Ethernet connection! Only the control element of the optional insulation monitor ELG may be connected here. Misuse can lead to damage to the isolating transformer as well as to the network!


Recycling of electrical and electronic equipment
Directive 2012/19/EU of the European Parliament and of the Council dated 04 July 2012 regarding waste electrical and electronic equipment

Symbols located on the packaging (if delivered as a spare part):

	Manufacturer of the isolating transformer
	Store dry
	Fragile! Handle with care
	Keep away from sunlight
	Temperature limits for the storage and transport of the isolating transformer in the original packaging
	Humidity range to which the isolating transformer may be exposed during storage and transport in the original packaging
	Air pressure range to which the isolating transformer may be exposed during storage and transport in the original packaging

Control panel of the ELG insulation monitor (optional)

1.4 Intended use

The intended use of the isolating transformer described in this manual is as follows:

Central power supply for medical electrical devices and devices tested to their relevant standards in patient environments and/or rooms used for medical purposes.
 Intended for fixed installation in the ITD "pro-cart" mobile equipment carriers

Forming a safe galvanic isolation of the supply and application circuit by use of a toroidal transformer with double and reinforced insulation between the primary and the secondary side

Limitation of the contact current to $\leq 100 \mu\text{A}$ in the normal state and $\leq 500 \mu\text{A}$ in single fault condition

Limitation of the earth leakage current to $\leq 5 \text{ mA}$ in the normal state and $\leq 10 \text{ mA}$ in single fault condition

Compliance with creepage and air distances as required according to standards

Compliance with EMC regulations and applicable standards

When medical electrical and in accordance with applicable standards-tested devices are combined to form a medical electrical system (MES) in the patient environment and/or in medically used rooms, the isolating transformer is connected upstream in accordance with the Regulation for Medical Devices (EU) 2017/745 to comply with the essential safety and performance requirements in accordance with Annex I MDR in line with IEC/EN 60601-1 (e.g: EDP accessories such as PC, monitor, printer, plotter, interfaces with medical electrical equipment as well as video equipment, diagnostic and monitoring cameras, diagnostic, measuring and testing devices and combinations thereof, etc.).

It should also be pointed out that the protective separation produced by the isolating transformer can be inadmissibly bridged by other electrical connections, such as computer network or USB connections to EDP devices arranged outside of the medical protection area. It is important to ensure that such connections are also equipped with an insulation location conforming to IEC/EN 60601-1 (e.g. Noratel LAN insulator MLI 1000).

1.5 Instructions for use

Prior to initial operation and installing the isolating transformer, please check the housing, power supply cord and mains plug for external injury or damage. If you find any injuries or damage, the isolating transformer (including the accessories) may not be put into operation under any circumstances. The repair or service is performed exclusively by the manufacturer or an authorised dealer of the mobile equipment carrier. The contact information can be found at the end of the operating manual.

Before you connect the isolating transformer to the supply mains, make sure that it is switched off. End devices may only be connected to the output of the isolating transformer if both the isolating transformer and the end device are switched off.

Note that the isolating transformer may only be operated on a supply mains with a functional protective conductor connection. The proper functioning of the protective conductor has to be ensured. In case of doubt, the supply mains has to be checked by a qualified specialist.

If the isolating transformer is operated in a room with potential equalisation, the isolating transformer has to be connected to the potential equalisation by means of a suitable cable. This is the only way to ensure that the user is not exposed to any dangerous touch currents.

1.6 Operating manual

The isolating transformer is equipped with an integrated, illuminated main switch. In addition, an external power switch, which is optionally located in the rail of the mobile equipment carrier, can be connected in series via an interface. If this is the case, the main switch located in the isolating transformer is permanently on and protected against accidental actuation by means of a cover hood. The mobile equipment carrier is then switched on and off exclusively via the external mains switch.

The input and output voltages set on delivery can be seen on the voltage selector switches (see section 1.0). To change the voltages, disconnect the isolating transformer from the supply mains and remove the plexiglass cover. Input and output voltages can be switched independently of each other. Replacement of the primary fuse inserts is only necessary if the supply voltage is changed. The corresponding values can be found in the table in section 1.9.

Make sure that the connected devices can be operated with the output voltage which is set at the isolating transformer. In case of doubt, please refer to the operating manual of the device manufacturer. Make sure that all end devices are off before they are connected to the isolating transformer.

Now, connect the desired end devices to the output sockets of the isolating transformer; then, follow the operating instructions of the mobile equipment carrier. If the isolating transformer is switched on, the set voltage is present at the output with a small time delay. The mains switch located in the isolating transformer, respectively in the rail lights up green. The end devices which are connected to the isolating transformer can now be switched on. Please observe the rules, regulations and requirements for the power supply cords, jumper cords and interface cables as well as their secure and tight fit.

The output power of the isolating transformer corresponds to the sum of the individual powers of the end devices which are connected to the isolating transformer and which are simultaneously in use. It is important to ensure that the required total power of the end devices does not exceed the maximum output power indicated on the isolating transformer nameplate. The available performance classes can be seen in Section 1.0.

1.7 Product features and designs

The following features apply to all available models and versions

- Protection Class IP 20 (in installed condition)
- Dimensions (L x W₁ x W₂ x H) 274 x 350 x 169 x 94 mm (660 und 1200 VA models)
- Dimensions (L x W₁ x W₂ x H) 274 x 350 x 169 x 120 mm (1600 VA models)
- POAG Bolt pursuant to DIN 42801 for the connection of the device to the potential equalisation of the building installation
- Toroidal transformer with a separate wired thermal switch 110 °C to trigger the overtemperature indicator circuit of the built-in earth leakage guard (only for devices of the ELG series)
- Input side electronic inrush current limiter with fast half-wave failure detection and extremely short duration of the network interruption of max. 50 milliseconds
- All-pole protection primary with fuse inserts
- Illuminated all-pole mains switch
- Primary and secondary separately switchable
- Plug protection against removal on the input side via clamping brackets
- Protection against removal on the output side via a corresponding mounting device on the ground frame of the mobile equipment carrier
- Protected against short-circuit and overload
- Classification as accessory for a medical device according to MDR (EU) 2017/745
- Class I medical device according to Rule 1 of Regulation (EU) 2017/745 according to Annex VIII, Chapter III, Paragraph 4.1
- Conformity and CE marking in accordance with Regulation (EU) 2017/745 as per Annex IX for Class I medical devices
- Conformity and approval according to EN 61558-1, EN 61558-2-4, IEC/ EN 60601-1, IEC/ EN 60601-1-2

1.8 Power supply cords, jumper cords

All power supply cords and jumper cords have to comply with the relevant standards and regulations of the individual countries in which the isolating transformers are used, e.g. UL/CSA/VDE/SEMKO/CHAR. H05VV-F3G1.5 with a maximum length of 5 meters is recommended as power supply cord. In the USA and Canada, a special mains connection cable type "SJ" or "SV" or better is required for use in hospitals. Connect the power supply cord to the isolating transformer on the appliance inlet and secure it with the protection against removal located on the transformer. Depending on the design of the plug, it may be necessary to exchange the fastening screws. Do not use screws larger than 50 mm in length; this could lead to damages inside the unit.

1.9 Changing the fuses

The fuse inserts in the input circuit protect the isolating transformer from short circuiting. If a fuse insert other than the one specified is used, there is a danger to the person or patient who is connected to the isolating transformer via the end device; moreover, serious damage to the connected end devices can also occur. The fuse inserts have to be marked with UL/CSA approvals for the American market, as well as with VDE/EN approvals for the European market. Please use time-lag fuses only (T label) for the isolating transformer.

As initial equipment Fuse inserts 5 x 20 mm pursuant to IEC 60127-2/5, UL 248-14, resp. CSA C22.2 no. 248.14 are used.

Please ensure that the isolating transformer is switched off before changing the fuse inserts and that no consumers are connected. Make sure the power cord is not connected to the power circuit. Then, use a screwdriver to open the fuse holders with a left turn. Where these are located is shown in Section 1.0. Insert the fuse inserts specified for the desired voltage and close the fuse holders by turning the cap clockwise. Reconnect the power cord and connect the isolating transformer to the consumers. First, switch on the isolating transformer and then the connected consumers.

Use only fuse inserts with the values listed below. Deviating values can impair the function or possibly lead to the destruction of the isolating transformer. Incorrect fuse values can also endanger the operating personnel and the patient.

The following fuse values are obligatory (time-lag, breaking capacity "H"):

Item No.	Drawing No.	Description	PRI 115 V	PRI 230 V
9-059-130128	ZV.9386.999	pro-cart 660 VA	6.30 A	3.15 A
9-059-130209	ZV.9575.999	pro-cart 660 VA ELG	6.30 A	3.15 A
9-059-130129	ZV.9387.999	pro-cart 1200 VA	12.50 A	6.30 A
9-059-130210	ZV.9576.999	pro-cart 1200 VA ELG	12.50 A	6.30 A
9-059-130212	ZV.9581.999	pro-cart 1600 VA	16.00 A	8.00 A
9-059-130213	ZV.9582.999	pro-cart 1600 VA ELG	16.00 A	8.00 A

DO NOT USE ANY FUSE INSERTS WITH HIGHER CURRENT VALUES!

1.10 Operating conditions

The following conditions have to be met for the operation of the isolating transformer:

- Installation in the designated equipment cart
- Maximum operating altitude above sea level 3000 m
- Ambient temperature range 10 °C to 40 °C ($t_a = 40$ °C)
- Humidity during operation 30 % to 75 %
- Air pressure during operation 700 hPa to 1060 hPa

1.11 Transport and storage conditions

The transport and storage conditions are stated on the packaging label of the isolating transformer. The required information is also listed here, in case the packaging can no longer be found:

- Store dry
- Fragile! Handle with care
- Keep away from sunlight
- Storage temperature limits -25 °C to + 70 °C
- Storage humidity range 5 % to 95 %
- Storage air pressure range 500 hPa to 1200 hPa

1.12 Cleaning, disinfection

The isolation transformer may not be dipped in liquids or autoclaved. No solvent-based chemical cleaners may be used. Before cleaning, the isolating transformer has to be disconnected from the mains supply.

The isolating transformer is not intended for operation in sterile rooms. It may not be submitted to any sterilisation procedures.

Tests have been carried out using the following products as examples:

Product	Manufacturer
Bacillol Plus	Bode
Cleanisept Wipes	Dr Schumacher
Mikrobac Tissues	Bode
Microzid Sensitive Wipes	Schülke
Terralin Protect	Schülke
Incidin Plus	Ecolab
Incidin Foam	Ecolab

Do not use any sharp-edged objects for cleaning. The isolating transformer may be put back into operation as soon as the cleaning agent has evaporated.

2.0 Description of the integrated insulation monitor (earth-leakage guard) ELG

In the case of the isolating transformers with an integrated isolation monitor, the evaluation electronic is installed in the housing, and the control and display panel is located in one of the bottom storage shelves. Both components are connected to each other by means of an interface cable laid in the vertical profile.

The description of the operating elements can be found in section 1.3.

2.1 Intended use

The ELG insulation monitor is used to monitor the insulation resistance of devices or device groups which are connected to the isolating transformer for protective separation. At the same time, the isolating transformer is monitored for its temperature behaviour. The evaluation is processor controlled.

The ELG operates in accordance with the specifications documented in DIN VDE 0107.

- Internal resistance > 100 k Ω
- Switching threshold > 50 k Ω
- Green operation – LED "On"
- Yellow signal – LED for signalling the limit value exceeded for the insulation resistance
- Yellow signal – LED for signalling of overtemperature
- Acoustic signalling of both
- Test buttons for function control of both functions
- Button for error acknowledgment

In addition, the following functions have been integrated:

- Function self-test of the ELG insulation monitor every time the device is switched on and in operation cyclically every 8 hours
- Error priority (acoustic) on insulation error

2.2 Product characteristics

- Modular design consisting of evaluation electronic in the transformer housing as well as operating and display panel in the bottom storage shelf
- Classification as accessory for a medical device according to MDR (EU) 2017/745
- Class I medical device according to Rule 1 of Regulation (EU) 2017/745 in accordance with Annex VIII, Chapter III, Paragraph 4.1
- Conformity and CE marking in accordance with Regulation (EU) 2017/745 as per Annex IX for Class I medical devices and VDE 0107
- Conformity pursuant to IEC/EN 60601-1, IEC/EN 60601-1-2, IEC/EN 60601-1-14
- Connection via CAT6 cable with RJ-45 plug connectors (pre-assembled and installed)
- Cable length approx. 3 m

2.3 Operating manual earth-leakage guard

Switch the main switch of the isolating transformer to "ON", the self-test of the ELG insulation monitor will start running automatically in the background within 5 seconds.

After completion of the self-test, the ELG insulation monitor is ready for operation, and the green LED is permanently on. The test is performed cyclically every 8 hours during operation and can also be triggered manually using the "Test" button.

The following test routine is performed during the manual test:

An insulation fault is simulated, the ISO LED is continuously lit, and a 2.4 kHz warning tone sounds continuously; both will extinguish after ca. 5 seconds.

Then, a temperature error is simulated, the TEMP LED is continuously lit and a 2.4 kHz warning tone sounds intermittently; both will extinguish after ca. 5 seconds.

An error can be detected as follows:

If an **ISOLATION ERROR** is detected, the ISO LED will light continuously, and a 2.4 kHz **warning tone** will sound **continuously**. The warning tone can be reset with the alarm-acknowledging button; the LED remains lit until the fault has been repaired.

If the isolation transformer is switched off and the fault is not repaired in the meantime, the above procedure starts again from the beginning.

When an insulation fault occurs, the acoustic and visual alarm remains on until they have been acknowledged:

First acknowledging: acoustic alarm off

Second acknowledging: visual alarm off

If a **TEMPERATURE ERROR** occurs, the TEMP LED lights continuously, and a 2.4 kHz **warning tone** sounds **intermittently**. The warning tone can be reset with the alarm acknowledging button; the LED remains lit until the fault has been repaired.

If the isolation transformer is switched off and the fault is not repaired in the meantime, the above procedure starts again from the beginning.

In the event of a simultaneous fault of insulation fault and overtemperature fault, the insulation fault always has priority with the acoustic alarm.

Self-test of the ELG insulation monitor:

In addition to the self-test which can be triggered by means of the "test" button, the ELG insulation monitor also performs a cyclic self-test in the period of ca. 8 hours; the test is also performed each time the unit is switched on.

The self-test lasts about 5 seconds and is not noticeable externally.

In the event of a fault, the green operating LED flashes with a frequency of 0.5 Hz, and the acoustic alarm sounds with the same frequency. The error messages cannot be reset with the acknowledging button.

2.4 Technical data of the ELG insulation monitor

Power supply	Nominal voltage	115/230 V ± 15 % (50/60 Hz)
Internal consumption		approx. 3 Watt
Response values (factory setting):	Response value	53 kΩ ± 2.5 kΩ
	Attention! The set value changes with the input voltage and is also dependent on its wave shape. The set value can deviate by up to ± 5%.	
Measuring circuit	Response time	> 0.5 s
	Test resistance	33 kΩ
	Internal resistance	
	AC	ca. 220 kΩ
	DC	ca. 320 kΩ
	Measuring voltage	12 V ± 0.3 V
	Max. measuring current	< 0.5 mA
Indicator elements	Max. external DC	250 V
	LED Green	Operation
	LED Yellow	Insulation fault
	LED Yellow	Temperature fault
Performed tests	Piezo	(2.4 kHz)
	Basic standard	DIN EN 55011: 2018-5
	Electrostatic Discharge	EN 61000-4-2: 2009-12
	Electromagnetic RF-Fields	EN 61000-4-3: 2011-4
	Fast Transients	EN 61000-4-4: 2013-1
	Surge	EN 61000-4-5: 2019-3
	RF conducted disturbance	EN 61000-4-6: 2014-8
	Voltage dips and interrupts	EN 61000-4-11: 2019-8

3.0 Elimination of faults

Do not attempt to repair the unit by yourself. If an improper repair attempt is determined, the guarantee will no longer be valid. For safety reasons, repairs and maintenance are performed exclusively by the manufacturer of the equipment cart.

No function:

1. Ensure that the ELG insulation monitor is connected to the isolating transformer, that this is connected to the mains, and that the ON-/ OFF switch is set to "ON".
2. Disconnect the entire system from the mains and remove all end devices from the output. Now check the fuse inserts of the isolating transformer.
3. Use or try another isolating transformer with an identical interface.
4. Check the circuit breaker of the supply circuit.
5. If you have any other problems, please contact the authorised dealer or manufacturer of the mobile equipment carrier.

In the case of mechanical or other damage, please immediately contact the manufacturer of the equipment cart directly:

ITD GmbH
Jahnstraße 1
DE 84347 Pfarrkirchen
Tel.: +49 89 614425-0
Fax: +49 89 614425-200
E-Mail: sales@itd-cart.com
www.itd-cart.com

Important note to users!

All serious incidents relating to the isolating transformer must be reported to the manufacturer and the regulatory authority of the member state in which the user is established.

3.1 Information about electromagnetic compatibility (EMC)

These isolating transformers have been evaluated pursuant to the following standards:

DIN EN 60601-1-2:2022-01 (IEC 60601-1-2:2014 + A1:2020)
 DIN EN 55011:2018-5
 DIN EN 61000-6-2:2019-11
 DIN EN IEC 61000-3-2:2023-10 n.a. acc.pt.7
 DIN EN 61000-3-3:2023-02
 DIN EN 61000-4-2:2009-12
 DIN EN IEC 61000-4-3:2021-11
 DIN EN 61000-4-4:2013-4
 DIN EN 61000-4-5:2019-3
 DIN EN 61000-4-6:2014-8
 DIN EN 61000-4-8:2010-11
 DIN EN IEC 61000-4-11:2021-10

Test conditions

Temperature	21.1 °C
Rel. humidity	51 %
Air pressure	930 – 1060 hPa
Power supply	230 V, 50 Hz

Evaluation and overview of the test results:

Interference Emissions		Result
Conducted Emissions	EN 55011, Group 1, Class B	Passed
Field strength	EN 55011, Group 1, Class B	Passed
Harmonics	EN 61000-3-2, Class A	Passed
Flicker	EN 61000-3-3	Passed
Interference immunity		
Electrostatic discharge	EN 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields	EN 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
Proximity fields from RF Wireless communication equipment	EN 61000-4-3	According Tab. 9 of IEC 60601-1-2:2014
Electrical fast transients (bursts)	EN 61000-4-4	± 2 kV 100 kHz repetition frequency
Surges	EN 61000-4-5	± 0,5 kV, ± 1 kV Line to line ± 0,5 kV, ± 1 kV, ± 2 kV line to ground
Conducted disturbances induced by RF fields	EN 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM-bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
RATED power frequency magnetic fields	EN 61000-4-8	30 A/m 50 Hz or 60 Hz
Voltage dips and interrupts	EN 61000-4-11	0 % UT; 0,5 cycle, 0 % UT; 1 cycle, 70 % UT; 25/30 cycles, 0 % UT; 250/300 cycle

3.2 Guarantee conditions

The guarantee for this isolating transformer lasts 24 months.

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© Sämtliche Inhalte und Texte sind urheberrechtlich geschützt. Nachdruck und Vervielfältigung nur mit ausdrücklicher Genehmigung durch den Hersteller.