INSTRUCTIONS FOR USE

Ambu® aScope™ Duodeno 2

For use by specialized healthcare professionals and experienced gastrointestinal endoscopists trained on ERCP procedures.


For use in professional healthcare facilities.

For use with Ambu® aBox™ 2.

Ambu

QUICK GUIDE

CONTENT

PAGE

English (Instructions for use)	04-19
Česky (Návod k použití)	20-36
Dansk (Brugsanvisning)	37-52
Deutsch (Bedienungsanleitung)	53-71
Ελληνικά (δηγίες Χρήσεως)	72-90
Español (Manual de instrucciones)	91-108
Suomi (Käyttöohje)	109-125
Français (Mode d'emploi)	126-143
Magyar (Használati útmutató)	144-161
Italiano (Manuale d'uso)	162-179
日本語 (使用法)	180-195
Nederlands (Gebruiksaanwijzing)	196-213
Norsk (Brukerveiledning)	214-230
Polski (Instrukcja obsługi)	231-248
Português (Manual de instruções)	249-266
Slovenčina (Návod na použitie)	267-284
Svenska (Instruktionshandbok)	285-300
Türkçe (Kullanım talimatları)	301-317
使用说明	318-331

1. Important information – Read before use

Read these instructions for use (IFU) carefully before using the aScope™ Duodeno 2. These instructions describe the function, setup and precautions related to operating the aScope™ Duodeno 2. Please be aware that these instructions do not explain or discuss clinical procedures.

Before initial use of the aScope™ Duodeno 2, it is essential for operators to have received sufficient training in clinical endoscopic techniques and to be familiar with the intended use, indications, warnings, cautions, and contraindications mentioned in these instructions.

There is no warranty on the aScope™ Duodeno 2.

Follow any valid official standards concerning the applicability of endoscopy and endoscopic treatment that have been established by hospital administration or other official institutions, such as academic societies on endoscopy. Before starting endoscopy and endoscopic treatment, thoroughly evaluate its properties, purposes, effects, and possible risks (their nature, extent, and probability). Even after starting endoscopy and endoscopic treatment, continue to evaluate its potential benefits and risks; stop endoscopy/treatment immediately and take appropriate measures if the risks to the patient ever outweigh the potential benefits.

The images must not be used as an independent diagnostic of any pathology. Physicians must interpret and substantiate any finding by other means and in the light of the patient's clinical characteristics. The instructions for use may be updated without further notice. Copies of the current version are available upon request.

In this document, “aScope™ Duodeno 2” refers to instructions that apply to the duodenoscope only and “aScope™ Duodeno 2 system” refers to information relevant to the aScope™ Duodeno 2, aBox™ 2 and accessories.

To enhance readability the aScope™ Duodeno 2 is referred to as endoscope or device and the aBox™ 2 is referred to as displaying unit within this IFU.

1.1. Intended use/Indication for use

The aScope™ Duodeno 2 is designed to be used with the aBox™ 2, endoscopic accessories (e.g., biopsy forceps) and other ancillary equipment (e.g., medical grade video monitor) for endoscopy and endoscopic surgery within the duodenum.

1.1.1. Intended patient population

Patients requiring endoscopy and/or endoscopic surgery within the duodenum.

1.1.2. Intended use environment

The device is intended to be used in professional healthcare facilities.

1.2. Intended user

For use by specialized healthcare professionals and experienced gastrointestinal endoscopists trained on ERCP procedures.

1.3. Contraindications

Do not use in patients with a bodyweight less than 15 kg.

1.4. Clinical benefits

Sterile, single use duodenoscope eliminates the risk of endoscope related cross contamination.

1.5. Warnings and cautions

WARNINGS

1. For single use only. Do not reuse, reprocess or resterilize, as these processes may leave harmful residues or cause endoscope malfunction. Reuse of the endoscope may cause cross-contamination potentially leading to infections.
2. Confirm that the opening of the insufflation/rinsing valve is not blocked or covered and that the insufflation pressure does not exceed the given limit. If gas is excessively insufflated into the patient, this may result in patient pain, bleeding, perforation and/or gas embolism.
3. Prior to use, always perform an inspection and functionality check according to sections 3.1 and 3.2. Do not use the endoscope if inspection or functionality check fails, as this can lead to patient injury or infection.
4. Patient leakage currents may be additive, when using energized endoscopic accessories. Do not use energized endoscopic accessories that are not classified as "type CF" or "type BF" applied part according to IEC 60601, as that could lead to patient leakage current that is too high.
5. Do not perform procedures with High Frequency (HF) endoscopic accessories if flammable or explosive gases are present in the gastrointestinal tract, as this may result in serious injury to the patient.
6. Always observe the live endoscopic image when inserting, withdrawing or operating the endoscope. Failure to do so may result in patient injury, bleeding and/or perforation. In the event of a temporarily impaired view or image loss, pause the procedure until view has been restored.
7. The temperature of the distal tip of the endoscope may exceed 41 °C (106 °F) due to endoscopic illumination. Avoid long periods of contact between the distal tip of the endoscope and the mucosa, as sustained contact may cause tissue damage.
8. Insertion or withdrawal of the endoscope while an endoscopic accessory is protruding from the distal end of the working channel may result in injury to the patient. Do not activate energized endoscopic accessories before they have been extended to an appropriate distance from the distal tip of the endoscope and can be seen in the endoscopic image. Otherwise, serious patient injury and/or endoscope damage may occur.
9. If the biopsy valve is damaged and/or left uncapped, it can reduce the efficiency of the endoscope's suction function and may pose a risk of infection due to leaking or spraying patient debris or fluids. When the valve is uncapped, cover it with a piece of sterile gauze to prevent leakage.
10. Always use gauze when pulling endoscopic accessories through the biopsy valve, as spraying or leaking patient debris or fluids could pose a risk of infection.
11. During the procedure always wear personal protective equipment (PPE) to protect against contact with potentially infectious material. Failure to do so may cause contamination potentially leading to infections.
12. Using HF endoscopic accessories with the endoscope may disturb the image on the displaying unit, which may lead to patient injury. To reduce disturbances, try alternative settings of the HF generator with lower Peak Voltage.
13. Portable Radio Frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the endoscope and the displaying unit, including cables specified by the manufacturer. Otherwise, degradation of the performance of the system could occur, which may lead to patient injury.
14. Only use the system with medical electrical equipment that complies with IEC 60601-1, and any applicable collateral and/or particular standards. Failure to do so may compromise medical system safety or lead to equipment damage.
15. Always lower the elevator before inserting or retracting the endoscope. Otherwise, the protruding elevator could cause patient injury.
16. Always use a patient mouthpiece in order to prevent the patient from accidentally biting the insertion tube, which could cause dental injury and/or damage to the endoscope and potentially dislodge parts/particles.

CAUTIONS

1. Handle with care: Do not drop, bump, bend, twist or pull any portion of the endoscope with excessive force as the endoscope may be damaged leading to failure in functionality.
2. Do not use excessive force to advance or retract endoscopic accessory (e.g., stents) or other material through the working channel. When inserting endoscopic accessories, confirm that the distal end of the endoscopic accessory is closed or completely retracted into its sheath. Failure to do so may cause damage to the endoscope.
3. Do not apply oil-based lubrication in the working channel as this may increase friction when inserting endoscopic accessories.
4. Do not coil the insertion tube or umbilical cord to a diameter of less than 12 cm (4.7 inches) as this may damage the endoscope.

1.6. Potential adverse events

Potential adverse events in relation to the Ambu® aScope™ Duodeno 2 (not exhaustive): Infection/inflammation (including but not limited to post-ERCP pancreatitis (PEP) and/or cholangitis), bleeding, tissue damage, perforation, thermal injuries, cardiopulmonary adverse events, air embolism, nausea, sore throat, abdominal pain and discomfort.


1.7. General notes

If, during the use of this device or as a result of its use, a serious incident occurs, please report it to the manufacturer and to your national authority.

2. Device description

The endoscope is to be connected to an Ambu displaying unit. For information about Ambu displaying units, please refer to the Ambu displaying units' instructions for use.

2.1. Device parts

aScope™ Duodeno 2 – Single-use device	Part numbers
	486001000 486004000


aScope™ Duodeno 2 might not be available in all countries. Please contact your local sales office.

Product name	Outer diameter	Inner diameter
aScope™ Duodeno 2	13.7 mm / 41.2 Fr	4.2 mm / 12.6 Fr

2.2. Product compatibility

The endoscope can be used in conjunction with:

- Ambu aBox™ 2.
- Insufflator for endoscopic gastrointestinal procedures with a constant flow of medical grade gas with a maximum supply pressure of 80 kPa (12 psi).
- Standard insufflation/rinsing management tubing sets compatible with Olympus endoscopes including sterile water bottle.
- A new disposable or sterilized insufflation/rinsing fluid management tubing set and sterile water bottle is required for each new procedure.
- Vacuum source to provide aspiration with a maximum vacuum of -76 kPa / -11 psi.
- Standard flexible suction tube.
- Irrespective of the chosen fluid management system, the suction container assembly used must feature overflow protection in order to prevent fluids from entering the system; this feature is commonly referred to as "self-sealing", "shut-off-filter" or similar.

No.	Part	Function
	<p>Control section</p> 	<p>The user holds the endoscope at the control section with the left hand.</p> <p>The control wheels and buttons may be manipulated with the left and/or right hand.</p>
1	Elevator control lever	<p>Turning this lever in the “U” direction raises the elevator shovel.</p> <p>Turning the lever in the opposite direction lowers the elevator shovel.</p>
2	Remote switches/ programmable buttons	<p>Switches/buttons to activate functions of the displaying unit, e.g., record video, image capture.</p> <p>The functions of the buttons, including the long and short press, are pre-configured from factory and can be reconfigured according to user’s preference.</p> <p>For detailed information, refer to the IFU of the displaying unit.</p>
3	Suction valve	<p>The valve is pressed down to activate suction to aspirate any fluids, debris or gas from the patient.</p> <p>The valve can be removed for unclogging or replacement.</p>
4	Insufflation/rinsing valve	<p>To activate insufflation, place a finger on the opening of the valve.</p> <p>To activate rinsing, press the valve fully down.</p>
5	Up/down control wheel	<p>When this wheel is turned in the “U” direction, the bending section moves up; when the wheel is turned in the “D” direction, the bending section moves down.</p>
6	Right/left control wheel	<p>When this wheel is turned in the “R” direction, the bending section moves right; when the wheel is turned in the “L” direction the bending section moves left.</p>
7	Up/down angulation lock	<p>Turning this lock in the “F” direction frees angulation.</p> <p>Turning the lock in the opposite direction locks the bending section at any desired position along the up/down axis.</p>
8	Right/left angulation lock	<p>Turning this lock in the “F” direction frees angulation.</p> <p>Turning the lock in the opposite direction locks the bending section at any desired position along the right/left axis.</p>
9	Biopsy valve	<p>Seals the working channel port.</p>

No.	Part	Function
10	Working channel port	The working channel port provides access to the working channel, which functions as: <ul style="list-style-type: none"> – Suction channel. – Channel for the insertion or connection of endoscopic accessories. – Fluid feed channel (from a syringe via the biopsy valve).
11	Insertion limit mark	Indicates the maximum length to which the endoscope may be inserted into the patient's body.
12	Insertion tube	The flexible insertion tube is inserted into the patient's upper gastrointestinal tract. The insertion tube has markings indicating the insertion depth. These markings do not represent an exact measurement of the insertion depth and serve only as an orientation.
13	Bending section	The bending section is the maneuverable part of the endoscope, that can be controlled by the control wheels and angulation locks.
14	Distal tip	The distal tip contains the camera, the light source, the working channel outlet and the insufflation/rinsing nozzle.
15	Camera	Enables visualization of the upper gastrointestinal tract and provides a live image stream.
16	Light source (LED)	Enables illumination of the upper gastrointestinal tract.
17	Elevator shovel	Manipulates endoscopic accessories and locks the guidewire.
18	Insufflation/rinsing nozzle	Nozzle for lens rinsing and insufflation.
19	Working channel outlet	This is the opening of the working channel at the distal end.
20	Suction connector	Connects the endoscope to the suction tube of the suction pump.
21	Insufflation/rinsing connector	Connects the endoscope to the sterile water bottle to enable insufflation and lens rinsing.
22	Lot number	Printed lot number of the endoscope.
23	Endoscope connector	Connects the endoscope to the grey connector port of the displaying unit. Ancillary equipment for suction, insufflation, lens rinsing can be attached to the endoscope connector.
24	Release button	Press the button when disconnecting the endoscope from the displaying unit.
25	Umbilical cord	Connects the control section with the endoscope connector.
26	Spare suction valve	Can be used to replace the existing suction valve in case of blockage or damage.

3. Use of aScope™ Duodeno 2

The number in grey circles refer to the Quick Guide on page 2.

Before each procedure, prepare and inspect a new endoscope as instructed below. Inspect other equipment to be used with this endoscope as instructed in their respective instruction manuals.

Should any irregularity be observed after inspection, follow the instructions as described in section 5, "Troubleshooting".

If the endoscope malfunctions, do not use it.

Contact your Ambu sales representative for further assistance.

3.1. Preparation and inspection of aScope™ Duodeno 2

Inspection of the device

Check that the pouch seal is intact, and that the device's expiration date has not yet been exceeded. In the case the pouch seal has been damaged or the expiration date has been exceeded, the endoscope must be discarded. **1a**

Carefully open the peel pouch and remove both the endoscope and the spare suction valve. **1b**

Make sure to use hygienic handling processes when completing the next steps.

Inspect the distal tip of the endoscope's insertion section for scratches, cracks, or other irregularities. **1e**

Turn the up-down and left-right control wheels in each direction until they stop and return to the neutral position. Confirm that the bending section functions smoothly and correctly, including that maximum angulation can be achieved and returned to the neutral position. **1c**

Operate the elevator control lever to confirm that the elevator shovel can be operated smoothly and correctly, and that it can be returned to its lowered position. **1d**

Confirm that the angulation locks are functional by locking and releasing the angulation locks. Turn the control wheels fully in all directions, lock angulation in a fully angulated position, and confirm that the bending section can be stabilized. Release the angulation locks and confirm that the bending section straightens out.

If applicable, confirm compatibility with applicable accessory devices.

A new endoscope should be readily available so the procedure can be continued in case a malfunction occurs.

Preparation for use

Prepare and inspect all necessary ancillary equipment as described in their respective IFU.

Power up the displaying unit. **2**

Carefully align the arrows on the connector of the endoscope cable with the grey port of the displaying unit to prevent damage to the connector.

Connect the endoscope to the displaying unit by plugging the endoscope connector into the corresponding grey port on the displaying unit. **3**

Check that the endoscope is securely locked to the displaying unit.

Attaching ancillary equipment

It is the responsibility of the user to follow all third-party manufacturer instructions and guidance applicable to the endoscopic fluid management and medical suction system chosen for use with the endoscope.

To perform patient examinations or procedures, maintain a safe working environment by ensuring that all fluid containers (e.g., sterile water bottle) are placed correctly and secured to prevent spillage.

Connecting insufflation/rinsing

Connect the endoscope using a disposable, high-level disinfected or sterilized insufflation/rinsing management tubing set. **4**

Please note that a new sterile or high-level disinfected water bottle must be used for each new procedure.

Check that the water level in the water bottle and remaining volume in the suction container are sufficient for the impending procedure.

Confirm that the fluid container's connection adapter fits properly and it cannot be rotated.

Connecting suction

As the endoscope is equipped with a standard suction connector, standard suction tubing can be used/is compatible. However, it is a requirement that the tubing and connector have a firm and tight connection.

When all other connections have been made, fit the end of the suction tubing securely over the suction connector located on the endoscope connector. **5**

Checking the image **6**

Verify that a live video image appears on the screen by pointing the distal tip of the endoscope toward an object, e.g., the palm of your hand.

Make sure that the orientation of the image is correct by moving the elevator control lever. The elevator shovel should move along the right-hand side of the video image.

If necessary, adjust image settings through the displaying unit settings and information menu.

If the image is fuzzy and/or blurry, wipe the lens using a moist sterile cloth.

Checking the remote switches/programmable buttons

Press every remote switch/programmable button and confirm that the specified functions work as expected.

Each remote switch/programmable button can be programmed to be sensitive on both short and long activation. See displaying unit IFU for further details.

Checking suction, rinsing and insufflation functionality

Press the suction valve fully down and confirm that the suction function works properly.

Cover the opening of the insufflation/rinsing valve and confirm that the insufflation function works properly.

Press the insufflation/rinsing valve fully down and confirm that the rinsing function works properly.

3.2. Operating the aScope™ Duodeno 2

Insertion of the endoscope **7**

Insert a suitable mouthpiece and place it between the patient's teeth or gums. If necessary, apply a medical-grade lubricant as indicated in section 2.2 to the insertion portion of the endoscope.

Before inserting the endoscope, make sure the elevator shovel is lowered before you insert the endoscope, and the bending section is straight.

Insert the distal end of the endoscope through the opening of the mouthpiece, then from the mouth to the pharynx while viewing the endoscopic image.

Do not insert beyond the proximal end maximum length mark.

Holding and maneuvering the endoscope

The control section of the endoscope is designed to be held in the operator's left hand.

The suction and insufflation/rinsing valves can be operated using the left index and middle fingers.

The up/down control wheel can be operated using the left thumb and supporting fingers.

The operator's right hand is free to manipulate the distal end via the insertion tube of the endoscope.

The right hand is intended to adjust the right/left control wheel and the angulation locks.

Angulation of the distal end

Operate the angulation control wheels as necessary to guide the distal end during insertion and observation.

The endoscope's angulation locks are used to hold the angulated distal end in position.

Insufflation/rinsing

Cover the opening of the insufflation/rinsing valve to feed gas from the insufflation/rinsing nozzle at the distal tip.

Press the insufflation/rinsing valve fully down to feed sterile water onto the objective lens.

Injection of fluids

Fluids can be injected through the working channel by inserting a syringe into the working channel port of the endoscope.

Insert the syringe completely into the port and press the plunger to inject fluid. Make sure you do not apply suction during this process, as this will redirect the injected fluids into the suction system.

Suction

Press the suction valve to aspirate excess fluids or other debris obscuring the endoscopic image.

For optimal suction capability it is recommended to remove endoscopic accessories entirely during suction.

Should the suction valve on the endoscope clog, remove and clean it or replace it with the spare suction valve attached on the mounting card.

Insertion of endoscopic accessories

Always make sure to select the correct size of endoscopic accessories for use in combination with the endoscope by consulting respective IFUs.

Accessories should be compatible if they are designed for working channels with an inner diameter (ID) of 4.2 mm / 12.6 Fr or less.

However, there is no guarantee that accessories selected using only this minimum working channel size will be compatible with the endoscope. Thus, compatibility of selected accessories should be assessed prior to the procedure.

Inspect the endoscopic accessory before use. Replace it if there is any irregularity in its operation or external appearance.

Confirm that the tip of the endoscopic accessory is closed or retracted into its sheath. Make sure to raise the elevator shovel by moving the elevator control lever in the "◀U" direction.

Hold the accessory approximately 4 cm / 1.5" from the biopsy valve and advance it slowly and straight towards the biopsy valve.

Insert the endoscopic accessory through the biopsy valve into the working channel using short strokes while observing the endoscopic image.

Confirm that the tip of the endoscopic accessory is touching the elevator shovel. Move the elevator control lever in the opposite direction to lower the elevator shovel. Advance the endoscopic accessory slightly and move the elevator control lever in the "◀U" direction.

Confirm that the endoscopic accessory is visible in the endoscopic image.

Manipulate the elevator control lever to adjust the height of the elevator shovel.

Withdrawal of endoscopic accessories

Ensure that the accessory is in a neutral position.

While gradually lowering the elevator shovel, withdraw the endoscopic accessory from the endoscope slowly through the biopsy valve.

If the accessory cannot be removed, retract the endoscope as described in the next paragraph while observing the endoscopic image.

Withdrawal of the endoscope 7

Stop using the image magnification (zoom) function of the displaying unit.

Move the elevator control lever in the opposite direction of the "◀U" direction until it stops.

Aspirate accumulated gas, blood, mucus or other debris by activating suction.

Move the up/down angulation lock to the "F" direction to release the angulation.

Turn the left/right angulation lock to the "F" direction to release the angulation.

Carefully withdraw the endoscope while observing the endoscopic image.

Remove the mouthpiece from the patient's mouth.

3.3. After use

Detach all tubes and tubing sets from the endoscope connector. 8

Press the disconnection button and disconnect the endoscope from the displaying unit. 9

Check the endoscope for any missing parts, evidence of damage, cuts, holes, sagging, or other irregularities on the insertion tube and bending section including the distal tip.

Should any irregularities exist, immediately determine if any parts are missing and take the necessary corrective action(s).

Discard the endoscope including all packaging and the spare suction valve in accordance with local guidelines for medical waste with electronic components. 10

Returning devices to Ambu

Should it be necessary to return an endoscope to Ambu for evaluation, please contact your representative at Ambu for instructions and/or guidance.

To prevent infection, it is strictly forbidden to ship contaminated medical devices.

As a medical device, the endoscope must be decontaminated on site prior to shipment to Ambu.

Ambu reserves the right to return contaminated medical devices to the sender.

4. Technical product specifications

4.1. aScope™ Duodeno 2 specifications

Insertion section dimensions	
Bending angle	Up: 120° Down: 90° Left: 90° Right: 110°
Max. insertion portion width [mm, ("), (Fr)]	13.7 / 0.5 / 41.2
Distal tip diameter [mm, ("), (Fr)]	13.7 / 0.5 / 41.2
Insertion tube outer diameter [mm, ("), (Fr)]	11.5 / 0.45 / 34.5
Working length [mm, (")]	1240 / 48.8
Min. working channel width [mm, ("), (Fr)]	4.2 / 0.17 / 12.6

Optics	
Field of view [°]	140
Direction of view [°]	6 (backward side viewing)
Depth of field [mm, (")]	5 – 100 / 0.19 – 3.93
Illumination method	LED
Maximum connection pressures	
Medical grade gas insufflator [kPa, (psi)]	Max. 80 / 12 (relative pressure)
Vacuum source [kPa, (psi)]	Max. -76 / -11 (relative pressure)
Operating conditions	
Temperature [°C, (°F)]	10 – 40 / 50 – 104
Relative humidity [%]	30 – 85
Atmospheric pressure [kPa, (psi)]	80 – 106 / 11 – 15
Sterilization	
Method of sterilization	Ethylene oxide (EO)
Storage and transportation conditions	
Transportation temperature [°C, (°F)]	-10 – 55 / 14 – 131
Storage temperature [°C, (°F)]	10 – 25 / 50 – 77
Transportation relative humidity [%]	10 – 95
Storage relative humidity [%]	10 – 85
Atmospheric pressure [kPa, (psi)]	50 – 106 / 7 – 15
Biocompatibility	
aScope™ Duodeno 2 is biocompatible	

5. Troubleshooting

The following tables show the possible causes of and countermeasures against challenges that may occur due to equipment setting errors or damage to the device.

Please contact your local Ambu representative if more detailed information is needed.

Prior to use, please do the pre-check as described in section 3.

Angulation and angulation lock

Possible problem	Possible cause	Recommended action
Increased resistance during control wheel operation.	The angulation lock is activated.	Release the angulation lock.
One or more of the control wheels do not turn.	Control wheel angulation locks are activated.	Release the angulation lock.
Angulation lock is not working.	Angulation lock is not correctly activated.	Activate the locking function by turning the angulation lock to the end stop.
Bending section does not angulate when control wheel is operated.	Endoscope is defective.	Withdraw the endoscope and connect a new one.

Possible problem	Possible cause	Recommended action
Max. bending angles cannot be reached.	Endoscope is defective.	Withdraw the endoscope and connect a new one.
Bending section angulates in the opposite direction.	Endoscope is defective.	Withdraw the endoscope and connect a new one.
Elevator shovel cannot be manipulated, or its movement is restricted.	Endoscope is defective.	Withdraw the endoscope and connect a new one.
Endoscopic accessory protrudes from the distal tip and cannot be withdrawn.	Elevator shovel raised.	Lower elevator shovel.

Rinsing/insufflation

Possible problem	Possible cause	Recommended action
Rinsing impaired or not possible.	Insufflation/rinsing fluid management tubing set not properly connected.	Connect the rinsing tubing properly to the endoscope.
	The water bottle is empty.	Replace the water bottle.
	Insufflator is not working, not switched ON or pressure setting is too low.	Refer to the insufflator IFU.
	Sterile water source setup insufficient.	Confirm that water source is installed according to its IFU.
	Insufflation/rinsing valve not fully activated.	Press the insufflation/rinsing valve fully down.
	Endoscope is defective.	Withdraw the endoscope and connect a new one.
Insufflation not possible or insufficient.	Insufflator is not connected, not switched ON or otherwise not working correctly.	Connect or switch ON compatible regulator. Adjust regulator settings. Refer to the insufflator IFU.
	Insufflation/rinsing fluid management tubing set not properly connected.	Connect the insufflation/rinsing fluid management tubing set properly to the endoscope.
	Sterile water source setup insufficient.	Confirm that water source is installed according to its IFU.
	Gas source is empty or remaining pressure too weak.	Connect a new gas source.
	Suction is activated.	Deactivate suction.
	Endoscope is defective.	Withdraw the endoscope and connect a new one.
Continuous insufflation without operating insufflation/rinsing valve.	Insufflation/rinsing valve opening is blocked.	Withdraw the endoscope and connect a new one.

Suction

Possible problem	Possible cause	Recommended action
Diminished or no suction.	Vacuum source/suction pump is not connected or not switched ON.	Connect vacuum source/suction pump and switch ON.
	Suction container is full or not connected.	Change the suction container if it is full. Connect a suction container.
	Suction valve is blocked.	Remove the valve and rinse with sterile water using a syringe and reuse the valve. Or replace the part with the spare suction valve.
	Biopsy valve is not properly connected.	Attach valve correctly.
	Biopsy valve cap is open.	Close cap.
	Vacuum source/suction pump too weak.	Increase vacuum pressure.
	Vacuum source/suction pump is defective.	Replace with a new vacuum source/suction pump.
	Working channel is blocked.	Flush sterile water with a syringe through the working channel.
	Endoscopic accessory inserted.	Remove endoscopic accessory.
	Endoscope is defective.	Withdraw the endoscope and connect a new one.
Continuous suction.	Suction valve remains pressed down.	Gently pull the suction valve up until suction stops.
	Suction valve defective.	Replace with the spare suction valve.

Working channel and use of endoscopic accessories

Possible problem	Possible cause	Recommended action
Working channel access is constricted or blocked (endoscopic accessories do not pass through channel smoothly).	Endoscopic accessory is not compatible.	Select a compatible endoscopic accessory.
	Endoscopic accessory is open.	Close the endoscopic accessory or retract it into its sheath.
	Working channel is blocked.	Try to unblock it by flushing sterile water into the working channel with a syringe.
	Biopsy valve is not open.	Open the cap of the biopsy valve.
	Incompatible lubricant used.	Withdraw the endoscope and connect a new one.
Insertion or withdrawal of endoscopic accessories is difficult.	High flexion of bending section.	Straighten the bending section as much as possible without losing the position of the endoscopic image.

Image quality and brightness



Possible problem	Possible cause	Recommended action
No video image.	Displaying unit or ancillary equipment is not switched ON.	Switch displaying unit and ancillary equipment ON.
	Endoscope connector is not properly connected to the displaying unit.	Connect the endoscope connector properly to the displaying unit.
	Endoscope is defective.	Withdraw the endoscope and connect a new one.
	Displaying unit is defective.	Contact your Ambu representative.
Image suddenly darkens.	Camera or illumination failure.	Switch ON LED as described in displaying unit IFU.
		Withdraw the endoscope and connect a new one.
Blurry image.	Objective lens is dirty.	Rinse the objective lens.
	Water drops on the outside of the lens.	Insufflate and/or rinse to remove water drops from the lens.
	Condensation on the inside of the lens during/after rinsing.	Pause rinsing until condensation disappears. Increase the water temperature in the water bottle and continue to use the endoscope.
	Displaying unit image settings improper.	See displaying unit IFU.
Flickering or disturbed images.	Signal interference from activated HF endoscopic accessory.	Use alternative mode or settings on the HF-generator with lower peak voltage (pV).
	Portable Radio Frequency (RF) communication equipment used too close to the device.	Move RF communication equipment further away or switch off.
	Signal interference due to fluoroscopy.	Stopping fluoroscopy will resolve any minor image interference it causes.
Dark or over-illuminated image.	Displaying unit image settings improper.	See displaying unit IFU.
	Endoscope is defective.	Withdraw the endoscope and connect a new one.
The color tone of the endoscopic image is unusual.	Endoscope is defective.	Withdraw the endoscope and connect a new one.
	Displaying unit image settings improper.	See displaying unit IFU.


Possible problem	Possible cause	Recommended action
Picture is frozen.	Endoscope is defective.	Withdraw the endoscope and connect a new one.
	Displaying unit is defective.	Contact your Ambu representative.
Unusual image contrast levels.	ARC mode unintentionally ON/OFF.	See displaying unit IFU.
	Improper image/ ARC settings.	See displaying unit IFU.

Remote switches/programmable buttons

Possible problem	Possible cause	Recommended action
The remote switches/programmable buttons are not working or not working properly.	The endoscope connector is not properly connected to the displaying unit.	Connect the endoscope connector properly to the displaying unit.
	Configuration of remote switches/programmable buttons changed by (other) user.	Return to standard configuration of the remote switches/programmable buttons or change the settings.
	Remote switch/programmable button is defective.	If feasible, operate the functions via the displaying unit (see displaying unit IFU); otherwise replace the endoscope.
	Displaying unit is defective.	Complete the procedure if possible; contact your Ambu representative.

6. Explanation of symbols used

Symbols for the aScope™ Duodeno 2 devices	Description	Symbols for the aScope™ Duodeno 2 devices	Description
	Working length of the insertion tube		Atmospheric pressure limitation
	Maximum insertion portion width (maximum outer diameter)		Humidity limitation
	Minimum instrument channel width (minimum inner diameter)		Temperature limit
	Manufacturer		Caution

Symbols for the aScope™ Duodeno 2 devices	Description	Symbols for the aScope™ Duodeno 2 devices	Description
	Field of view		Country of manufacture. Date of manufacture
	Do not use if package is damaged		Global trade item number
	Do not re-use		UL mark on electronic products (UL Recognized Component Mark for Canada and the United States)
	Consult instructions for use		Use-by date
	Batch Code		Catalogue number
pMax CO₂	Maximum relative supply pressure by insufflator. Values are depicted in kPa/psi		Medical Device
	Single sterile barrier system	pMax VAC	Maximum relative negative pressure supplied by vacuum source. Values are depicted in kPa/psi
	Indicates that the product is in compliance with European legislation for medical devices and has passed conformity assessment by EU Notified Body		UK Conformity Assessed. Indicates that the product is in compliance with UK legislation for medical devices and has passed conformity assessment by UK Approved Body
	UK Responsible Person		Importer

A full list of symbol explanations can also be found on ambu.com/symbol-explanation.

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Ambu and other trademarks are trademarks of Ambu A/S.
Ambu A/S
Ambu Ltd