

CARDIOLOGY ACCESSORIES

Features and Benefits

- Fixes electrode leadwires on patients during Holter and Stress test
- Reduces movement artefacts on the ECG
- Minimizes the risk of electrodes falling off
- MR safe
- Adhesive size: Ø 55 mm
- PVC free
- Latex free
- Biocompatible



Skin Fix Ref: 2703/50

Units/pouch	Units/inner box	Units/outer box
50	500	3000

Features and Benefits

- Allows skin preparation before electrode placement
- Gently abrasive to reduce skin impedance
- Adhesive on one side, to use with one finger
- Small size fits finger tip
- Size: 19 x 28 mm
- PVC free



Skin Prep Ref: 2121M

Units/pouch	Units/inner box	Units/outer box
500	1500	15000

CARDIOLOGY ACCESSORIES

Features and Benefits

- Connects 3 mm pin cables to tab electrodes
- Ensures secure connection
- Easy to use
- Reusable
- PVC free



Snap Clip Ref: 0334M

Units/pouch	Units/inner box	Units/outer box
10	N/A	1000

Features and Benefits

- Connects 4 mm banana pin to snap or tab electrodes
- Ensures secure connection
- Allows quick and easy one-hand application
- Reusable
- PVC free



Multi Clip Ref: 2952/10

Units/pouch	Units/inner box	Units/outer box
10	N/A	N/A



Hambu A/S Baltorpbakken 13 2750 Ballerup Denmark T +45 7225 20 00 F +45 7225 2053 ambu.com