Ambu® WhiteSensor WSP30

ECG Electrode - Single use



Key Benefits

- Small/medium-sized electrode for pediatric as well as adult use
- Highly conductive solid gel and strong adhesion for a stable and secure tracing
- Flexible foam backing

Ambu® WhiteSensor WSP30

The small/medium-sized Ambu WhiteSensor WSP30 features a highly conductive solid gel with strong adhesion to ensure a good signal quality during pediatric and adult ECG monitoring.

Thanks to the flexible foam backing material and the teardrop shape, the electrode ensures a good fit on most patients, ease of use and comfort during application.





Solid gel



Specifications

Dimensions				
Electrode size (W x L in mm)	35 x 30 511 1			
Adhesive area (in mm²)				
Height excluding connector/ wire (in mm)				
Sensor				
Sensor material	Silver/silver chloride (Ag/AgCl)			
Gel system	Solid gel			
Sensor area (in mm²)	44			
Gel area/measuring area (in mm²)	201			
Electrical data (ANSI/AAMI)				
AC impedance	≤ 150 Ω			
DC offset voltagel	≤ 1.0 mV			
Defibrillation overload recoveryl	≤ 15 mV			
Rate of change of polarization potential	≤ 0.3 mV/s			
Impedance after SDR (ACZ)	≤ 150 Ω			
Combined offset instability and internal noise	≤ 90 µV			
Bias current tolerance (over 8 hours)	≤ 15 mV			

Environment		
PVC-free electrode	Yes	
Electrode not made with natural rubber latex	Yes	
PVC-free packaging	Yes	
X-Ray & MRI		
Radiolucent	No	
MR Conditional	No	
Shelf life		
Opened pouches	1 month	
Unopened pouches	24 months*	

^{*} From date of production.

Materials

8 hours)

Electrode		
Bio-compatible	Yes	
Sensor	Silver/silver chloride (Ag/AgCI)	
Outer carrier	Polyethylene foam (enhanced)	
Outer adhesive	Polyacrylate	
Connector	Stainless steel stud	
Liner	PET	
Packaging		
Pouch (outer/middle/ inner layers)	Polyester/Aluminum/ Polythene	
Boxes	Corrugated cardboard	

Options		
Connectors	S	
Precaution		
Single use only		

Available configurations

Reference number	Units/card	Units/pouch	Units/inner box	Units/outer box
WSP30-00-S/50	5	50	600	7200

Product

WSP30-00-S • Wire length

Connector

Distributor in EU Ambu A/S

Baltorpbakken 13 DK - 2750 Ballerup Denmark Tel. +45 7225 2000 ambu.com

Manufacturer **MEDICO ELECTRODES** INTERNATIONAL LTD. [AN ISO13485 COMPANY] Plot 142A/11, 12, 27, 28 & 29, Noida Special Economic Zone, Noida, Uttar Pradesh - 201305, India www.medicoelectrodes.com