DATASHEET

The Ambu WhiteSensor 40713 features a skinfriendly conductive solid gel and small size to ensure a good signal quality during ECG monitoring applications in paediatrics.

The electrode has breathable and soft cloth backing material.


It is radiolucent and MR Conditional.

Key benefits

- Skin-friendly solid gel supports the quality of the signal
- Soft cloth backing allows the skin to breathe
- Small size is ideal for paediatric patients
- Radiolucency allows the electrode to remain on the patient during X-ray procedures
- MR Conditional

Ambu[®] WhiteSensor 40713

ECG Electrode Single use





Solid gel



Cloth backing



Silver sensor



Radiolucent



MR Conditional



Breathable backing



SPECIFICATIONS

Dimensions		
Electrode size (length x width in mm)	38 x 19	
Skin contact size (length x width in mm)	38 x 19	
Adhesive area (in mm²)	722	
Sensor		
Sensor material	Silver/silver chloride (Ag/AgCl)	
Gel system	Solid gel	
Gel area (in mm²)	722	
Electrical data (ANSI/AAMI)		
AC impedance - typical	560 Ω	
DC offset potential - typical	4.6 mV	
Defibrillation overload recovery - typical	19.67mV	
Rate of change for polarisation potential - typical	-0.39 mV/s	
Combined offset instability and internal noise	3.3 μV	
Bias current tolerance	10.1 mV	

Environment	
Electrode is not made with natural	
rubber latex	

Packaging is PVC-free

X-ray & MRI		
Radiolucent	Yes	
MR Conditional	Yes*	
Shelf life		
Opened pouch	1 month	
Unopened pouch	24 months**	

*Static magnetic field of 1.5 Tesla and 3 Tesla only. Maximum spatial gradient field of 25.000 Gauss/cm or 250 Tesla/m. Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

** from date of production

CONNECTORS



S/RT = Graphite Stud

MATERIALS

Electrode	
Bio-compatible	Yes
Backing material	Non-woven (Cloth)
Connector	Graphite snap-type
Release liner	Siliconised polyester
Packaging	
Pouch (outer/middle/inner layers)	Polyester/Aluminium/Polyethylene
	,
Boxes	Cardboard
Boxes Options	, ,
····	, ,
Options	Cardboard

AVAILABLE CONFIGURATIONS*

Reference number	Units/liner	Units/pouch	Units/inner box	Units/outer box
40713	3	3	300	3000

*some configurations may not be available worldwide

Ambu

M Ambu A/S

Baltorpbakken 13 2750 Ballerup Denmark T +45 72 25 20 00 ambu.com

