Ambu[®] WhiteSensor 4530

ECG Electrode - Single use



Key Benefits

- Highly conductive wet gel
- Flexible foam backing
- 3 mm fitting connector
- Radiolucent
- MR Conditional
- High-tack adhesive foam
- Biocompatible

Ambu[®] WhiteSensor 4530

The Ambu® WhiteSensor 4530 features a highly conductive wet gel system supporting a fast and reliable trace during 12-lead ECG procedures. Presented conveniently with 5 pieces to a card. Thanks to the strong high-tack adhesive foam, the electrode supports workflow by minimizing electrode drop-off. The 3 mm slot connector which is radiolucent makes it possible for the patient to undergo X-ray and scanner procedures.









Radiolucent



MR Conditional



White Sensor

Specifications

Dimensions

Dimensions			
Electrode size (length x width in mm)	32 X 44		
Skin contact size (length x width in mm)	32 X 41		
Adhesive area (in mm²)	1079		
Height excluding connector (in mm)	1.5		
Sensor			
Sensor material	Silver/silver chloride (Ag/AgCl)		
Gel system	Wet gel		
Gel area (in mm²)	201		
Sensor area (in mm²)	77		
Electrical data (ANSI/AAMI)			
AC impedance - typical	200 Ω		
DC offset potential - typical	1.0 mV		
Defibrillation overload recovery - typical	12 mV		
V	0.3 mV/s		
Combined offset instability and internal noise	60 µV		
Bias current tolerance	5 mV		

Environmen

Environment	
PVC-free electrode	Yes
Electrode not made with natural rubber latex	Yes
PVC-free packaging	Yes
X-Ray & MRI	
Radiolucent	Yes
MR Conditional	Yes*
Shelf life	
Opened pouch	1 month
Sealed pouch	24 months**

* Static magnetic field of 1.5 Tesla and 3 Tesla only. Maximum spatial gradient field of 25.000 Gauss/cm or 250 Tesla/m. Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning.

** From date of production.

Materials

Electrode	
Bio-compatible	Yes
Coopera	Polyester/Polyurethane/
Sponge	Reticulated
Backing material	Polyethylene foam (PE)
Backing material adhesive	Polyacrylate
Supporting label	Polyethylene
Connector	3 mm graphite slot
Release liner	Siliconised polyester

Packaging	
Pouch (outer/middle/ inner layers)	Polyester/Aluminum/ Polythene
Boxes	Corrugated cardboard
Options	
Connectors	A
Precaution	
Single use only	
Accessories	
Skin Prep	Ref: 2121M

Available configurations

4530 5 60	300	2400

Distributor in EU Ambu A/S

Baltorpbakken 13 DK - 2750 Ballerup Denmark Tel. +45 7225 2000 Fax +45 7225 2053 www.ambu.com

Manufacturer MEDICO ELECTRODES INTERNATIONAL LTD.

[AN ISO13485 COMPANY] Plot 142A/11, 12, 27, 28 & 29, Noida Special Economic Zone, Noida, Uttar Pradesh - 201305, India www.medicoelectrodes.com

CE US: Rx only