# Ambu® WhiteSensor 4500M

# ECG Electrode - Single use



#### **Key Benefits**

- Highly conductive wet gel
- Good adhesion
- Flexible foam backing
- Radiolucent
- MR Conditional

### Ambu® WhiteSensor 4500M

The Ambu WhiteSensor 4500M features a wet gel with good adhesion to ensure a good signal quality during short-term ECG monitoring applications.

Thanks to the flexible foam backing material, the electrode ensures ease of use and comfort during application. It is radiolucent and MR Conditional.









MR Conditional





## Specifications

Dimensions				
Electrode size (length x width in mm)	50 x 48			
Adhesive area (in mm²)	1584			
Height excluding connector (in mm)	1			
Sensor	Sensor			
Sensor material	Silver/silver chloride (Ag/AgCI)			
Gel system	Wet gel			
Gel area (in mm²)	201 79			
Sensor area (in mm²)				
Electrical data (ANSI/AAMI)				
AC impedance - typical	95 Ω			
DC offset potential - typical	1.0 mV			
Defibrillation overload recovery - typical	12 mV			
Rate of change for polarisation potential - typical	0.3 mV/s			
Combined offset instability and internal noise	9 μV			
Bias current tolerance	5 mV			

Environment		
PVC-free electrode	Yes	
Electrode not made with natural rubber latex	Yes	
PVC-free packaging	Yes	
X-Ray & MRI		
Radiolucent	Yes	
MR Conditional	Yes*	
Shelf life		
Opened pouch	30 days	
Sealed pouch	24 months**	

<sup>\*</sup> Static magnetic field of 1.5 Tesla and 3 Tesla only.

Maximum spatial gradient field of 25.000 Gauss/cm or
250 Tesla/m. Maximum whole body averaged specific
absorption rate (SAR) of 2 W/kg for 15 minutes of scanning.

#### **Materials**

Electrode	ectrode		
Bio-compatible	Yes		
Sponge	Polyester Polyurethane		
	Reticulated		
Backing material	Polyethylene foam (PE)		
Backing material adhesive	Polyacrylate		
Supporting label	Polyethylene		
Release liner	Polyester (PET)		
Connector	Carbon filled ABS		

Packaging	kaging			
Pouch (outer/middle/ inner layers)	Polyester/Aluminum/ Polythene			
Boxes	Corrugated cardboard			
Options				
Connectors	S			
Precaution				
Single use only				
Accessories				
Skin Preps	Ref: 2121M			

## Available configurations

Reference number	Units/card	Units/pouch	Units/inner box	Units/outer box
4500M	3	30	300	3000

Distributor in EU Ambu A/S

Baltorpbakken 13 DK - 2750 Ballerup Denmark Tel. +45 7225 2000 Fax +45 7225 2053 www.ambu.com Manufacturer MEDICO ELECTRODES INTERNATIONAL LTD.

[AN ISO13485 COMPANY]
Plot 142A/11, 12, 27, 28 & 29,
Noida Special Economic Zone,
Noida, Uttar Pradesh - 201305, India
www.medicoelectrodes.com

<sup>\*\*</sup> From date of production.