Ambu[®] WhiteSensor 4831Q

ECG Electrode - Single use



Key Benefits

- Solid gel
- Flexible foam backing
- Ergonomic leaf shape
- Convenient presentation of 3 pieces/strip
- Small size can be used on children and adults

Ambu® WhiteSensor 4831Q

Adult and paediatric short term foam monitoring electrode with solid gel. The electrode is made in the leaf shape and is conveniently presented with 3 pcs to a strip. Thanks to the flexible foam backing material and the ergonomic shape, the electrode ensures ease of use and comfort during application.







Specifications

Dimensions			
Electrode size	37.5 x 37.5		
(length x width in mm)	57.5 × 57.5		
Skin contact size (mm²)	1306		
Adhesive area (in mm²)	992		
Height excluding connector (in mm)	1		
Sensor			
Sensor material	Silver/silver chloride		
	(Ag/AgCl)		
Gel system	Solid gel		
Gel area (in mm²)	201		
Sensor area (in mm²)	79		
Electrical data (ANSI/AAMI)			
AC impedance - typical	140 Ω		
DC offset voltage - typical	1.0 mV		
Defibrillation overload recovery - typical	11 mV		
Rate of change of polarization	0.2		
potential - typical	0.2 mV/s		
Combined offset instability	10.00		
and internal noise	28 µV		
Bias current tolerance	5 mV		

Environment			
PVC-free electrode	le Yes		
Electrode not made with natural rubber latex	Yes		
PVC-free packaging	Yes		
X-Ray & MRI			
Radiolucent	No		
MR Conditional	No		
Shelf life			
Opened pouches	1 month		
Unopened pouches*	24 months*		

*from date of production

Materials

Yes				
Silver/silver chloride (Ag/AgCl)				
Polyethylene foam (PE)				
Polyacrylate				
Stainless steel				
Polyethylene				
PET, printed				
Packaging				
Polyester/Aluminum/Polythene				
Corrugated cardboard				

S	
	S

Available configurations*

Reference number	Units/card	Units/pouch	Units/inner box	Units/outer box
4831Q	3	60	600	7200

*some configurations may not be available worldwide

Distributor in EU Ambu A/S

Baltorpbakken 13 DK - 2750 Ballerup Denmark Tel. +45 7225 2000 ambu.com Manufacturer MEDICO ELECTRODES INTERNATIONAL LTD.

[AN ISO13485 COMPANY] Plot 142A/11, 12, 27, 28 & 29, Noida Special Economic Zone, Noida, Uttar Pradesh - 201305, India www.medicoelectrodes.com

