Ambu® WhiteSensor 4841P

ECG Electrode - Single use



Key Benefits

- Solid gel
- · Flexible foam backing
- Ergonomic leaf shape
- Radiolucent
- MR Conditional

Ambu® WhiteSensor 4841P

Adult short term foam monitoring electrode with solid gel. The electrode is made in the leaf shape and is conveniently presented with 3 pcs to a strip. Thanks to

the flexible foam backing material and the ergonomic shape, the electrode ensures ease of use and comfort during application. It is radiolucent and MR Conditional.







Foam backing



Radiolucent



MR Conditional





Specifications

Dimensions		
Electrode size	42 x 42	
(length x width in mm)		
Skin contact size	1624	
(length x width in mm)		
Adhesive area (in mm²)	1302	
Height excluding connector (in mm)	1	
Sensor		
Sensor material	Silver/silver chloride	
	(Ag/AgCl)	
Gel system	Solid gel	
Gel area (in mm²)	201	
Sensor area (in mm²)	57	
Electrical data (ANSI/AAMI)		
AC impedance - typical	150 Ω	
DC offset voltage - typical	1.0 mV	
Defibrillation overload recovery - typical	11 mV	
Rate of change of polarization	0.3 mV/s	
potential - typical	0.5 1117/5	
Combined offset instability	28 μV	
and internal noise	20 μν	
Bias current tolerance	5 mV	

Environment				
PVC-free electrode	ree electrode Yes			
Electrode not made with natural rubber latex	Yes			
PVC-free packaging	Yes			
X-Ray & MRI				
Radiolucent	Yes			
MR Conditional	ditional Yes*			
Shelf life				
Opened pouches	1 month			
Unopened pouches	24 months**			

*Static magnetic field of 1.5 Tesla and 3 Tesla only. Maximum spatial gradient field of 25.000 Gauss/cm or 250 Tesla/m. Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning.

**from date of production

Materials

Electrode	
Bio-compatible	Yes
Sensor	Silver/silver chloride (Ag/AgCl)
Backing material	Polyethylene foam (PE)
Backing material adhesive	Polyacrylate
Connector	Carbon filled ABS
Supporting label	Polyethylene
Release liner	PET, printed
Packaging	
Pouch (outer/middle/inner layers)	Polyester/Aluminum/Polythene
Boxes	Corrugated cardboard

Option		
Connectors	S	
Precaution		
Single use only		

Available configurations*

Reference number	Units/card	Units/pouch	Units/inner box	Units/outer box
4841P	3	60	600	7200

*some configurations may not be available worldwide

Distributor in EU Ambu A/S Baltorpbakken 13 DK - 2750 Ballerup Denmark Tel. +45 7225 2000 ambu.com Manufacturer
MEDICO ELECTRODES
INTERNATIONAL LTD.
[AN ISO13485 COMPANY]
Plot 142A/11, 12, 27, 28 & 29,
Noida Special Economic Zone,
Noida, Uttar Pradesh - 201305, India
www.medicoelectrodes.com

CE US: Rx only