

DATASHEET

This Infant IV Arm has a special, extremely thin, synthetic skin, and rubber tubing with appropriately small lumen and thin walls.

The cephalic and basilic vein are accessible, as well as the dorsal venous arch on the hand.

This arm is the same as on the Complete/Deluxe Infant CRiSis™ Manikins.
Fluid supply stand not included.

Infant IV Arm



KEY BENEFITS

- Simulating a realistic infant arm
- Training of intravenous access
- Special, extremely thin, synthetic skin
- Tubing with small lumen and thin walls
- Allows access on cephalic and basilic vein, and dorsal venous arch on the hand

ORDERING SPECIFICATIONS

Item no.

LF03637U

Infant IV Arm

incl. infant arm, 2 IV bags with clamps, 1 pint blood, winged infusion set

Ambu

■ Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com