DATASHEET

An excellent choice for nasal endoscopy, flexible laryngoscopy and swallowing assessments.

The single-use aScope™ 4 RhinoLaryngo Slim video endoscope is designed with a small 3.0 mm outer diameter to help minimize patient discomfort during procedures. It delivers high quality imaging, is always ready and requires no reprocessing or repairs.

Excellent imaging helps you perform procedures confidently

The combination of the high image quality from the aScope™ 4 RhinoLaryngo Slim with the full-HD Ambu® aView™ 2 Advance monitor delivers excellent imaging. You can count on clear sharp images, which make it easy to identify anatomical structures.

Single-use efficiency streamlines your workflow

Whether you're on a consult or in the clinic, time is always of the essence. aScope 4 RhinoLaryngo is always available and easy to take with you streamlining your entire workflow.

Ambu® aScope™ 4 RhinoLaryngo Slim



KEY BENEFITS

High quality imaging on the full-HD monitor Ambu® aView™ 2 Advance

It's an always available solution that supports increasing patient volumes

Improves workflow with no hassle of cleaning and repair

Wide range of connectivity options

A brand-new endoscope with the same high quality at every procedure




SPECIFICATIONS

Optical system	
Field of view	85°
Direction of view	0° (forward viewing)
Depth of field	6 - 50 mm
Illumination method	LED
Insertion portion	
Bending section	130° up, 130° down
Insertion cord diameter	3.0 mm (0.12")
Distal tip diameter	3.5 mm (0.14")
Maximum diameter of insertion portion	3.5 mm (0.14")
Working length	300 mm (11.8")
Storage	
Temperature*	10 - 25 °C (50 to 77 °F)
Relative humidity	10 - 85 %
Atmospheric pressure	50 - 106 kPa
Transportation	
Temperature	-10 - 55 °C (14 to 131 °F)
Relative humidity	10 - 95 %
Atmospheric pressure	50 - 106 kPa
Operating environment	
Temperature	10 - 40 °C (50 to 104 °F)
Relative humidity	30 - 85 %
Atmospheric pressure	80 - 106 kPa
Altitude	≤ 2000 m
Sterilisation	
Method of sterilisation	ETO

ORDERING SPECIFICATIONS**

Item no.	Quantity
510001000 - aScope 4 RhinoLaryngo Slim	5





Ambu® aScope $^{\text{™}}$ 4 RhinoLaryngo Slim and Ambu® aView $^{\text{™}}$ 2 Advance



aScope $^{\mathsf{M}}$ 4 RhinoLaryngo Slim is always comfortable to hold thanks to its lightweight and ergonomic design



■ Ambu A/S

Baltorpbakken 13 2750 Ballerup Denmark +45 72 25 20 00 ambu.com