

Ambu A/S
 Baltorpbakken 13
 DK-2750 Ballerup
 Denmark
 T +45 72 25 20 00
 F +45 72 25 20 50
 www.ambu.com

Ambu is one of the world's leading suppliers of innovative diagnostic and life supporting devices, carefully designed to provide genuine value and ultimately improve the quality of patient care. Our largest business areas are Respiratory Care, Cardiology and Neurology in which the most important products are ventilation products for artificial respiration and single-use electrodes for ECG tests and neurophysiological mappings.

Headquartered in Copenhagen, Denmark, approx. 98% of Ambu's products is sold on export markets worldwide through Ambu's foreign sales companies or distributors. The total number of employees is 1216 of whom 313 work in Denmark and 903 in countries outside Denmark.

Ambu is listed on the Copenhagen Stock Exchange.

AMBU BUSINESS AREAS

- RESPIRATORY CARE
- CARDIOLOGY
- NEUROLOGY
- TRAINING
- IMMOBILIZATION



Ambu[®] AuraFlex



Always part of you anaesthesia team

THE NEED: SAFETY, FLEXIBILITY AND EASE OF WORK

The challenge of a flexible laryngeal mask is the combination of needs:

The ease of insertion, fast and accurate positioning and security of the airway needed by the anesthesiologist and the flexibility and room to work needed by the surgeon.

Also, the protection of the patients from residual proteins is crucial when using a flexible laryngeal mask, as this is often directly exposed to blood and other secretes due to its position within the surgical area.

THE SOLUTION: AMBU AURAFLEX

Ambu AuraFlex combines these needs in an innovative, single-use solution.

Ambu AuraFlex is so flexible that it can easily be taped well away from the surgical field, improving surgical access without any loss of seal. The flexible tube is also reinforced with kink-free wire, which eliminates the risk of airway tube occlusion.

Furthermore, the integration of pilot tube in the airway tube makes handling and insertion easy and dramatically reduces the risk of cutting the pilot tube during surgery.

Also, Ambu AuraFlex includes some of the same well-known features as the rest of the Aura family, such as the reinforced tip for fast and accurate positioning, extra soft, thin cuff with easy glide surface for a quick and smooth insertion and teeth marks for a visual check of the mask position.

And finally, the single-use concept eliminates the risk of cross-contamination and at the same time eliminates all costs and hassle of cleaning and sterilization.

AuraFlex - the new standard in disposable flexible laryngeal masks.



FEATURES AND BENEFITS AT A GLANCE

- Disposable, flexible laryngeal mask specially designed for ENT, ophthalmic, dental and other head and neck surgeries.
- High flexibility of the airway tube enables positioning away from the surgical field without loss of seal
- Integrated pilot tube eases the surgeons work and reduces the risk of cutting the pilot tube during surgery
- Kink-free wire-reinforced airway tube eliminates the risk of airway tube occlusion
- Single-use, eliminates the risk of cross contamination and the cost and hassle of cleaning and sterilization
- Full range of sizes from 2 through to 6 covering all needs
- Thin and soft feel of the cuff with EasyGlide surface for easy and smooth insertion and removal
- Reinforced tip reducing the risk of the tip folding during insertion, thereby ensuring easy, fast, and accurate positioning
- Depth marks for visual control of correct positioning of the mask
- Packaged sterile and ready for use. The cuff is protected during transport and storage by a specially designed disposable cuff protector
- Color-coded packaging with directions for use