

Ambu® ET Tube Holder



Key Benefits

- Designed to secure the ET tube in place after insertion in the trachea
- Comfortable closed cell foam padded backing
- Eliminates the need for tape and timely taping procedures
- Integrated bite block
- Quick Connect Ratchet System
- Non-latex strap with strap inducer for easy application around head
- Packaged 20 units per case
- Single patient use
- MR safe
- Latex free
- PVC free
- CT Scan safe

Ambu® ET Tube holder

Ambu's endotracheal tube holder allows the rescuer to apply holder with one hand in just seconds. The Ambu ET Tube Holder is easy to use and can be placed on patients needing a 3 mm - 9 mm ET Tube.

For use on both adults and children.

Materials

Product	Material
Blue strap	Silicone
White strap	Nylon
Inducer (yellow)	Polypropylen (PP)
Tube body holder	Polycarbonat
Packaging	
Pouch	Polyethylen (PE)
Box	Cardboard

Ordering specifications

Item no.	Description
320 264 040	Ambu ET Tube Holder (Blue Rubber Strap)
320 264 041	Ambu ET Tube Holder (White Soft Velcro Strap)

Packaging	
Quantity/pouch	1
Quantity/box	20

