

Ambu[®] aFlow

CO₂ sampling line




Ambu aFlow CO₂ sampling line maintain reliable results during sidestream CO₂ measurement. Available in different lengths of 2-3 meters.

Key Benefits

- Disposable, Single Patient Use
- Medical grade tubing
- 100 % Flow/occlusion tested
- Not made with natural rubber latex
- Gas sampling lines used to monitor CO₂, O₂, and NO levels

CO₂ sampling line



Ambu item no.	350500100	350500200	350500300
Quantity / Box	100	100	100
Material: PE / EVA / PVC	✓	✓	✓
Tube length	3 m	2 m	3 m
Tube size	1.2 x 2.8 mm	1.2 x 2.8 mm	1.2 x 2.8 mm
Newport connector	-	✓	
Connection mechanism	Male - Male	Male - Male	Male - Male

Distributed by:
Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com

Manufactured by:
**Altera Tibbi Malz.
 San. ve Tic. A.S.**
 Ibni Melek OSB Mh. Tosbi Yol 4 Sok. No:29
 35900 Tire Izmir
 TURKEY
 Tel: (+90 232) 513 51 10 (Pbx)
 www.altera.com.tr