


Ambu® aFlow

Catheter Mounts



Extendible Catheter Mounts



Ambu aFlow offers a wide range of catheter mounts covering several size and material specifications. Catheter Mounts are available in adults and pediatrics sizes and as smoothbore, and extendible in the following material; PVC and PP. Extendible catheter mounts provides a lightweight and flexible solution for clinicians to adjust the catheter position as they see fit. Ambu aFlow catheter mounts are suitable for connections between airway devices and breathing circuits and compatible with all types of breathing and ventilator circuits following relevant standards.


Key Benefits

- Suitable for connection between invasive airway devices and breathing circuits
- Smoothbore products provide low resistance to flow while external coil reduces kinks and blockages. Moreover, they prevent fluid building up and thereby reducing risks of bacterial growth thanks to the smooth inner surface.
- The extendible PP products provide a lightweight and flexible solution and allows clinician to adjust the position of the catheter mount
- Some catheter mounts are equipped with 360 degree double swivel elbow connector heads which increase ease of use
- Double cap helps to maintain PEEP during suctioning or bronchoscopy.
- Luer lock port for gas sampling
- 15M or 22F adapters available



Ambu item no.	350300300	350300500	350300600	350300700	350300800	350300900
Size	Adult / Pediatric	Adult / Pediatric	Adult / Pediatric	Adult / Pediatric	Adult / Pediatric	Adult / Pediatric
Quantity/Box	50	50	30	20	20	50
Material	PP	PP	PP	PP	PP	PP
Tube length	15 cm	15 cm	15 cm	15 cm	15 cm	15 cm
Tube I.D. size	15 mm	15 mm	15 mm	15 mm	15 mm	15 mm
Double Swivel Elbow	✓	✓	✓	✓	✓	✓
Ports	Double port	Double port	W/out ports	Double port	W/out ports	Double port
Red safety cap	✓	✓	✓	✓	✓	✓
Swivel Head with Soft Cap	✓	✓	-	✓	✓	✓
Bacterial viral filter	-	-	-	With soft cap	With soft cap	HME filter with soft cap
Connection patient end	15F	15F	15F	15F	15F	15F
Connection circuit end	15M	22F	22F	22 M	22M	22M
Compressible Volume ml	Closed: 25 ml, Open: 40 ml	Closed: 25 ml, Open: 40 ml	Closed: 25 ml, Open: 40 ml	93 ml	73 ml	Closed: 103 ml, Open: 117 ml

Catheter Mount Extension Tubes



Ambu item no.	350301300	350301400	350301500	350301600	350301700	350301800	350301900
Size	Adult / Pediatric	Adult / Pediatric	Adult	Adult	Adult	Adult	Adult
Quantity/Box	50	50	20	20	20	20	20
Tube Type	Smoothbore	Extendible	Extendible	Extendible	Extendible	Extendible	Extendible
Material	PVC	PP	PP	PP	PP	PP	PP
Tube length	15 cm	15 cm	15 cm	15 cm	15 cm	15 cm	15 cm
Tube I.D. size	15 mm	15 mm	15 mm	15 mm	15 mm	15 mm	15 mm
Red safety cap	✓	✓	✓	✓	✓	✓	✓
Bacterial viral filter	-	-	HME with soft cap	With soft cap	With soft cap	HME with soft cap	HME with soft cap
Connection patient end	15F	15F	15F	15F	15F	15F	15F
Connection circuit end	22F	22F	22M	22M	22M	22M	22M
Compressible Volume ml	27 ml	Closed: 22 ml, Open: 37 ml	76 ml	56 ml	56 ml	78 ml	78 ml

Distributed by:
Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com

Manufactured by:
Alterta Tibbi Malz.
San. ve Tic. A.S.
 İbni Melek OSB Mh. Tosbi Yol 4 Sok. No:29
 35900 Tire İzmir
 TURKEY
 Tel: (+90 232) 513 51 10 (Pbx)
 www.alterta.com.tr