

Ambu[®] aFlow

Adult Corrugated Anaesthesia Circuits




Corrugated anaesthesia circuits are lightweight and effective circuits. They have directly bonded connectors to the tube which helps in avoiding leakage risks.

Key Benefits


- Directly bonded connectors reducing leakage risks
- Numerous variants of kits with and without watertraps to fit every need
- The soft, lightweight and flexible materials enhance a secure seal while still allowing easy connection and disconnection

Adult Corrugated Anaesthesia Circuits



Ambu item no.	350200200	350200300	350200400	350200500
Size	Adult	Adult	Adult	Adult
Quantity/Box	20	15	20	20
Material	PE	PE	PE	PE
Tube length	1.6 m	2.4 m	1.6 m	1.6 m
Tube I.D.	22 mm	22 mm	22 mm	22 mm
Elbow Connector with luer cap	✓	✓	✓	✓
Fixed Elbow Connector Head	✓	✓	✓	✓
Y-connector Specifications	Dual ports	W/out ports	W/out ports	Swivel Y-connector with luer cap port
Red safety cap	✓	✓	✓	✓
Straight connector 22-22MM	✓	✓	-	-
Compliance @60 cmH ₂ O	1.69	2.959	1.698	1.72
Resistance to Flow @30 L/min	0.013	0.014	0.0125	0.023

Adult Corrugated Anaesthesia Circuits Kits



Ambu item no.	350200600	350200700	350200800	350200900	350201000	350201100	350201900	350202000
Size	Adult	Adult	Adult	Adult	Adult	Adult	Adult	Adult
Quantity/Box	15	10	15	15	15	15	15	15
Material	PE	PE	PE	PE	PE	PE	PE	PP
Tube length	1.6 m	1.8 m	1.8 m	2.4 m	2.4 m	3.6 m	1.6 m	1.6 m
Tube I.D.	22 mm	22 mm	22 mm	22 mm	22 mm	22 mm	22 mm	22 mm
Elbow Connector with luer cap	✓	✓	✓	✓	✓	✓	✓	✓
Fixed Elbow Connector Head	✓	✓	✓	✓	✓	✓	✓	✓
Y-connector Specifications	Dual port	W/out ports	W/out ports	W/out ports	W/out ports	W/out ports	Dual port	Dual port
Red safety cap	✓	✓	✓	✓	✓	✓	✓	✓
Bacterial viral filter	-	-	-	-	-	-	-	✓
Double water trap	-	-	-	-	-	-	✓	✓
Extention limb	0.8 m	0.8 m	1.8 m	0.8 m	1.2 m	1.2 m	0.8 m	1.2 m
Breathing bag, 2 L Volume	-	-	2 L	2 L	2 L	2 L	-	2 L
CO ₂ Sampling line	-	-	-	-	3 m	3 m	-	3 m
Tube Clip	-	-	-	-	-	-	✓	-
Compliance @60 cmH ₂ O	1.49	1.859	1.99	2.44	3.37	2.93	1.82	1.94
Resistance to Flow @30 L/min	0.0125	0.0135	0.0125	0.014	0.0135	0.0165	0.016	0.0165

Distributed by:
Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com

Manufactured by:
Alterta Tibbi Malz.
San. ve Tic. A.S.
 İbni Melek OSB Mh. Tosbi Yol 4 Sok. No:29
 35900 Tire Izmir
 TURKEY
 Tel: (+90 232) 513 51 10 (Pbx)
 www.altera.com.tr