


Ambu® aFlow

Extendible Anaesthesia Circuits




Extendible circuits may be extended up to 3 times their packed length, which allows the system to be set to the desired length and position, optimizing the compressible volume. Available in adult 22 mm and pediatric 15 mm sizes.

Key Benefits


- Maintain position set by healthcare professional minimizing strain on airway connections
- Tubes may be extended up to 3 times their packed length
- Tubes can be set to a desired length and position optimizing the compressible volume.

Pediatric Extendible Anaesthesia Circuits



Ambu item no.	350205300	350205400	350205500
Size	Pediatric	Pediatric	Pediatric
Quantity/Box	35	35	15
Material	PP	PP	PP
Tube length	1.8 m	1.5 m	3.0
Tube I.D.	15 mm	15 mm	15 mm
Elbow Connector with luer cap	✓	✓	✓
Fixed Elbow Connector Head	✓	✓	✓
Y Connector Specifications	W/out port	Dual port	W/out port
Detachable Y-connector	-	-	-
Red safety cap	✓	✓	✓
Extention limb		0.75 m	1.5 m
Breathing bag, Volume	-	-	1 L
Compliance @60 cmH ₂ O	0.92	0.75	1.25
Resistance to Flow @30 L/min	0.023	0.0205	0.0305

Adult Extendible Anaesthesia Circuits



Ambu item no.	350203000	350203100	350203300	350203400	350203500	350203600
Size	Adult	Adult	Adult	Adult	Adult	Adult
Quantity/Box	35	35	35	30	30	20
Material	PP	PP	PP	PP	PP	PP
Tube length	1.8 m	1.6 m	1.8 m	1.8 m	2.0 m	1.8 m
Tube I.D.	22 mm	22 mm	22 mm	22 mm	22 mm	22 mm
Elbow Connector with luer cap	✓	✓	✓	✓	✓	✓
Fixed Elbow Connector Head	✓	✓	✓	✓	✓	✓
Y Connector Specifications	W/out port	Dual port	W/out port	W/out port	W/out port	W/out port
Detachable Y-connector	-	✓	-	-	-	-
Red safety cap	✓	✓	✓	✓	✓	✓
Bacterial viral filter	-	-	-	-	-	✓
Water trap	-	-	-	-	-	Double
Extention limb	-	-	0.9 m	0.9 m	1.5 m	1.5 m
Breathing bag, Volume	-	-	-	-	2 L	2 L
CO ₂ Sampling line	-	-	-	3 m	-	3 m
Compliance @60 cmH ₂ O	1.25	1.17	1.25	1.22	1.43	1.32
Resistance to Flow @30 L/min	0.022	0.02	0.021	0.021	0.023	0.0265

Distributed by:
Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com

Manufactured by:
Altera Tibbi Malz.
San. ve Tic. A.S.
 İbni Melek OSB Mh. Tosbi Yol 4 Sok. No:29
 35900 Tire Izmir
 TURKEY
 Tel: (+90 232) 513 51 10 (Pbx)
 www.altera.com.tr

