

Ambu[®] aFlow

Adult Smoothbore Anaesthesia Circuit kits




Smoothbore anaesthesia circuits have low compliance and transparent tubing which makes viewing the inside of the tube easier. The inner smooth surface allows for easy drainage of condensed water and contributes to reduce the risk of any interruptions during anaesthetic procedures.

Key Benefits

- Lower tube compliance
- Transparent tubing which makes viewing inside easier
- Smooth inner surface for easy drainage of condensed water
- Soft flexible connectors ensure a perfect seal with the anaesthesia machine inlet and outlet.

Adult Smoothbore Anaesthesia Circuit kits



Ambu item no.	350202100	350202200	350202300	350202400	350202500	350202600	350202700	350202800	350202900
Size	Adult	Adult	Adult	Adult	Adult	Adult	Adult	Adult	Adult
Quantity/Box	20	15	15	15	10	15	15	15	10
Material	PVC	PVC	PVC	PVC	PVC	PVC	PVC	PVC	PVC
Tube length	1.8 m	1.8 m	1.8 m	1.8 m	1.8 m	1.8 m	1.8 m	1.8 m	1.8 m
Tube I.D.	19 mm	19 mm	19 mm	19 mm	19 mm	19 mm	19 mm	19 mm	19 mm
Elbow connector with luer cap	✓	✓	✓	✓	✓	✓	✓	✓	✓
Fixed elbow connector head	✓	✓	✓	✓	✓	✓	✓	✓	✓
Y-connector dual port	-	✓	✓	✓	✓	✓	✓	✓	✓
Detachable y-connector	-	✓	-	-	✓	-	-	-	-
Red safety cap	✓	✓	✓	✓	✓	✓	✓	✓	✓
Bacterial viral filter	-	-	-	✓	HME	✓	✓	✓	HME
Water trap	-	-	-	-	Single	Double	Double	Double	Double
Extention limb	-	1.2 m	1.2 m	1.2 m	0.6 m	0.6 m	1.2 m	1.2 m	0.6 m
Breathing bag, volume	-	-	2 L	2 L	-	-	2 L	2 L	-
CO ₂ Sampling line	-	-	-	-	-	-	3 m	2 x 3 m	-
Catheter mount	-	-	-	-	Double swivel, Extensible PP 15 mm	-	-	-	Double swivel, Extensible PP 15 mm
Compliance @60 cmH ₂ O	1.15	1.25	1.25	1.19	1.31	1.32	1.34	1.24	1.23
Resistance to flow @30 L/min	0.018	0.0135	0.0175	0.0165	0.02	0.02	0.189	0.02	0.02

Distributed by:
Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com

Manufactured by:
Altera Tibbi Malz.
San. ve Tic. A.S.
 İbni Melek OSB Mh. Tosbi Yol 4 Sok. No:29
 35900 Tire İzmir
 TURKEY
 Tel: (+90 232) 513 51 10 (Pbx)
 www.altera.com.tr