

# Ambu<sup>®</sup> Disposable Face Mask



## Key Benefits

- Crystal clear dome for easy observation of the patient's condition
- Dome with thumb rest for an easy grip
- Soft, shaped, airfilled cuff allows a tight fit to the patient's face
- Supplied with hook ring; hook ring can easily be removed if not needed
- The single patient use concept offers minimized risk of cross infection
- All sizes are delivered individually packed, non-sterile, in packs of 20 units
- Phthalate-free material

## Ambu<sup>®</sup> Disposable Face Mask

### Description

The Disposable Face Mask is an anatomically correct shaped face mask designed for oxygenating and ventilating any patient.

The face mask comes in a full range of 6 different sizes ensuring that the user can always find a mask that fits any patient.

### Features

It is a single patient use face mask that comes both with and without a check valve. The face mask has a thumb rest for easy handling while providing a tight seal.

The mask shape ensures that the inner shape of the dome is spacious enough for the mouth and nose also when pressure is applied to the mask.

## Ordering specifications

Item no.	Type
000 252 051	Ambu Disposable Face Mask Neonate
000 252 052	Ambu Disposable Face Mask Infant
000 252 053	Ambu Disposable Face Mask Toddler
000 252 054	Ambu Disposable Face Mask Small adult/Child
000 252 055	Ambu Disposable Face Mask Medium adult
000 252 056	Ambu Disposable Face Mask Large adult
000 252 951	Ambu Disposable Face Mask w. checkvalve Neonate
000 252 952	Ambu Disposable Face Mask w. checkvalve Infant
000 252 953	Ambu Disposable Face Mask w. checkvalve Toddler
000 252 954	Ambu Disposable Face Mask w. checkvalve Small adult/Child
000 252 955	Ambu Disposable Face Mask w. checkvalve Medium adult
000 252 956	Ambu Disposable Face Mask w. checkvalve Large adult



Ambu Disposable Facemasks with checkvalve.

## Materials

Product	Phthalate-free	Material
Dome	YES	Polyvinylchloride
Cuff	YES	Polyvinylchloride
Check valve	YES	Polycarbonate
Packaging	Quantity	Material
Pouch	1	PP (Polypropylene)
Box	20	Cardboard



Ambu Disposable Facemasks connected to Ambu SPUR II.

## Size

Type	Weight
Neonate	9 g
Infant	14 g
Toddler	25 g
Small adult/Child	33 g
Medium adult	47 g
Large adult	49 g

### Ambu A/S

Baltorpbakken 13  
 DK - 2750 Ballerup  
 Denmark  
 Tel. +45 7225 2000  
 Fax +45 7225 2053  
 www.ambu.com