

Ambu® Open Cuff Silicone Face Mask



Key Benefits

- Transparent for visual check of bleeding, vomitus, and spontaneous breathing
- Available in 6 different sizes
- All models can be autoclaved repeatedly at 134° C

Ambu® Open Cuff Silicone Face Mask

Ambu Open Cuff Silicone Face Masks are designed for use with manual and automatic resuscitators and ventilators.

The masks are made of soft silicone and shaped for a comfortable fit for both adults and children.

Ordering specifications

Item no.	Type
000 251 001	Ambu Open Cuff Silicone Face Mask size 0 infant
000 251 002	Ambu Open Cuff Silicone Face Mask size 1 Infant
000 251 003	Ambu Open Cuff Silicone Face Mask size 2 Child
000 251 004	Ambu Open Cuff Silicone Face Mask size 3 Child
000 251 005	Ambu Open Cuff Silicone Face Mask size 4 Adult/child
000 251 006	Ambu Open Cuff Silicone Face Mask size 5 Adult

Materials

Product	Material
Mask	Silicone rubber

Size

Type	Weight
0 Infant	13 g
1 Infant	26 g
2 Child	41 g
3 Child	54 g
4 Adult/Child	76 g
5 Adult	102 g



Ambu Open Cuff Silicone Face Mask sizes 0, 1 and 2



Ambu Open Cuff Silicone Face Mask sizes 3, 4 and 5