

Ambu[®] Silicone Face Mask



Key Benefits

- Transparent dome for visual check of bleeding, vomitus and spontaneous breathing
- Thumb rest for a tight, easy seal
- Comfortable design and texture
- The buttoned-on cuff enables easy removal for cleaning
- Can be autoclaved repeatedly at 134°C (except size oA)
- Available in 7 sizes
- Size o with boring impedes ventilation pressures higher than approx. 40 cm H₂O, enabling the user to administrate correct ventilation on small children with a body weight of at least 15 kg (3 years)

Ambu[®] Silicone Face Mask

The Ambu Silicone Face Mask is designed for use with manual and automatic resuscitators and ventilators. Ambu Silicone Face Mask is self-inflating for an easy and tight seal. The cuff pressure may be regulated (except for mask size oA).

Ordering specifications


Item no.	Type
000 313 000	Ambu Silicone Face Mask size oA infant
000 314 000	Ambu Silicone Face Mask size o Infant w.boring
000 315 000	Ambu Silicone Face Mask size o Infant
000 316 000	Ambu Silicone Face Mask size 2 Child
000 312 000	Ambu Silicone Face Mask size 3/4 Child
000 317 000	Ambu Silicone Face Mask size 5 Adult
000 319 000	Ambu Silicone Face Mask size 6 Adult

Materials

Product	Material
Cuff	Silicone rubber
Cuff size oA	Silicone rubber with foam filling
Dome	Polysulphone
Mask connector	Silicone rubber

Size

Type	Weight
oA Infant	28 g
o Infant w.boring	39 g
o Infant	38 g
2 Child	63 g
3/4 Adult/Child	78 g
4 Adult	80 g
5 Adult	105 g
6 Adult	105 g



Ambu Silicone Face Mask connected to Ambu Mark IV Resuscitator



Ambu Silicone Face Mask connected to Ambu Oval Plus Silicone resuscitators