

Ambu® UltraSeal



Key Benefits

- Very soft anatomically shaped cuff enabling a tight seal with minimum applied pressure
- Flexible, self-expanding dome
- Shoulder grip fitting different hand sizes
- Crystal clear dome for easy observation of the patient's condition
- Supplied with a clear hook ring hook ring can easily be removed if not needed
- The single patient use concept offers minimized risk of cross infection
- All sizes are delivered individually packed in a transparent, easy to open bag
- Color coding on mask label for fast and easy identification of size
- All sizes are delivered individually packed, non-sterile, in packs of 20 units
- Phtalate-free material

Ambu® UltraSeal

Ambu UltraSeal is disposable face mask designed specifically for anesthetic departments.

The anatomically shaped cuff secures a tight seal and enables optimal ventilation.

Complete line with 4 different sizes from size 2-3 Toddler to size 6 Large adult. All sizes come with and without check valve.

Ordering specifications

Item no.	Description
● 305 010 000	Ambu UltraSeal # 4 adult/Child
● 305 012 000	Ambu UltraSeal # 5 Medium adult
● 305 014 000	Ambu UltraSeal # 6 Large adult
● 305 107 000	Ambu UltraSeal # 3 Toddler, with check valve
● 305 011 000	Ambu UltraSeal # 4 Small adult/Child, with check valve
● 305 013 000	Ambu UltraSeal # 5 Medium adult, with check valve
● 305 015 000	Ambu UltraSeal # 6 Large adult, with check valve

Size

Type	Weight without check value	Weight with check value
# 3 Toddler	37 g	38 g
# 4 Small adult/Child	44 g	45 g
# 5 Medium adult	53 g	54 g
# 6 Large adult	72 g	73 g

Materials

Product	Phthalate-free	Material
Dome	YES	Polyvinylchloride
Cuff	YES	Polyvinylchloride
Check valve	YES	Polycarbonate
Packaging	Quantity	Material
Pouch	1	PP (Polypropylene)
Box	20	Cardboard



Ambu UltraSeal without check valve



Ambu UltraSeal with check valve