

Ambu® UltraSeal paediatric



Key Benefits

- **Transparent and clear dome**
Secures an easy observation of the patient's condition during use
- **Soft and flexible cuff**
The cuff is soft and provides good comfort while creating a tight seal
- **Single patient use**
No need for time-consuming cleaning, sterilizing and disinfecting procedures. The added risk of cross contamination is reduced
- **Low weight and good edge grip**
Optimal user handling and support during use
- **Phthalate-free material**

Ambu® UltraSeal paediatric

Description

The UltraSeal face mask is an anatomically correct shaped facemask designed specifically for anaesthesia departments.

The facemask comes in four different sizes ensuring that the user can always find a mask that fits any paediatric patient.

Features

It is a single patient use face mask that comes both with and without a check valve. The face mask has good edge grip which makes it easy to use while providing a tight seal. The mask shape ensures that the inner shape of the dome is spacious enough for the mouth and nose also when pressure is applied to the mask.

Ordering specifications

Item no.	Description	Unit pr. box
● 305 102 000	UltraSeal # 0	20
● 305 104 000	UltraSeal # 1	20
● 305 106 000	UltraSeal # 2	20
● 305 108 000	UltraSeal # 3	20
● 305 101 000	UltraSeal # 0 with checkvalve	20
● 305 103 000	UltraSeal # 1 with checkvalve	20
● 305 105 000	UltraSeal # 2 with checkvalve	20
● 305 107 000	UltraSeal # 3 with checkvalve	20



Ambu UltraSeal without check valve

Size

Type	Weight without check value	Weight with check value
# 0	9.1 g	9.9 g
# 1	12.4 g	13.2 g
# 2	19.2 g	20 g
# 3	25.3 g	26.1 g



Ambu UltraSeal with check valve

Materials

Product	Phthalate-free	Material
Dome	YES	Polyvinylchloride
Cuff	YES	Polyvinylchloride
Check valve	YES	Polycarbonate
Packaging	Quantity	Material
Pouch	1	PP (Polypropylene)
Box	20	Cardboard