

# Ambu® LifeKey



## Key Benefits

- Semi-transparent face shield and one-way valve allowing rescuer to visually check for patient's lip colour and vomitus
- Prevents contact with patient's saliva or blood
- Suitable for both adults and children
- Elastic ear bands fit over patient's ears to keep the mask in place
- Can be carried on a key chain, in a pocket, a small handbag and a first-aid kit

## Ambu® LifeKey

The Ambu LifeKey is a proven, easy to use, and inexpensive barrier with a one-way valve that prevents direct contact with patient's mouth, nose and face.

The Ambu LifeKey helps protect rescuer and increases bystanders' motivation for CPR in emergency situations.

Standard resuscitation techniques can be used without modification.

Comes with Ambu logo or as private labelling. Available in hardcase, softcase or plastic bag.

## Ordering specifications

Ambu LifeKey	
248 204 102	Holster Pack containing lifekey, one set of gloves and alcohol pad
248 210 001	Refill, set of 25 masks in a plastic bag
248 204 101	Lifekey in Hardcase with Ambu logo White
248 204 104	Lifekey in Hardcase with Ambu logo Red
248 206 102	Lifekey in Softcase with Ambu logo Black
248 206 103	Lifekey in Softcase with Ambu logo Blue
248 206 104	Lifekey in Softcase with Ambu logo Red
248 206 105	Lifekey in Softcase with Ambu logo Yellow
248 206 106	Lifekey in Softcase with Ambu logo Green
248 206 108	Lifekey in Softcase with Ambu logo Orange
Ambu LifeKey private label	
248 204 301	Lifekey in Hardcase with private label White
248 204 303	Lifekey in Hardcase with private label Blue
248 204 304	Lifekey in Hardcase with private label Red
248 204 305	Lifekey in Hardcase with private label Yellow
248 204 306	Lifekey in Hardcase with private label Green
248 206 303	Lifekey in Softcase with private label Blue
248 206 304	Lifekey in Softcase with private label Red
248 206 305	Lifekey in Softcase with private label Yellow
248 206 306	Lifekey in Softcase with private label Green
248 206 308	Lifekey in Softcase with private label Orange

## Materials


Product	Material
Mask	Polyethylene
Softcase	Nylon
Hardcase	Styrolux

## Size

Type	
Folded	3 x 2 x 1 cm
Hardcase	3.7 x 3.7 x 1.5 cm
Softcase	5 x 5 x 1.5 cm
Weight app	10 g
<b>Hardcase</b>	
Printing area (logo size)	20 x 20 mm Available in colours: black, blue, green, red, white and yellow
<b>Softcase</b>	
Printing area (logo size)	14 x 39 cm Available in colours: black, blue, red, yellow, green, orange and white
<b>Order</b>	
Minimum	250 pieces



Ambu LifeKey unfolded



Ambu LifeKey in use



Ambu LifeKey private label