

# Ambu® Perfit® Extrication Collars



## Key Benefits

- Ambu collars are fully disposable. You no longer have to worry about collars being returned to the source of origin and cross contamination issues are avoided
- Ambu collars come with a pre-shaped flip-chin piece to ensure they can be stored completely flat. Furthermore the flip-chin piece facilitates intubation and removal of vomit from a patient's mouth
- One-piece extrication collar simplifies application and at the same time there are no extra pieces that can get lost
- Nasal cannula holders are located on the collar body to secure the oxygen tubes (except for infant)
- All collars are radiotranslucent and come with a ventilated posterior shell for fluid drainage and a large hole for easy access. All of these features allow for patient examinations to take place safely without removing the collar

## Ambu® Perfit® Extrication Collars

The Ambu Perfit collars are designed to assist with the maintenance of neutral alignment, prevention of lateral sway and anterior-posterior flexion and extension of the spine during transport and routine patient care or movement.

## Specifications

Description	
Shell material	Polyethylene
Sizes	1 (infant), 2 (paediatric), 3 (neckless), 4 (short), 5 (regular) and 6 (tall)
Colour	White/light grey

## Ordering specifications

Item no.	Description
264 201 000	Perfit Infant: size 1 / cartons of 20 pcs.
264 202 000	Perfit Paediatric: size 2 / cartons of 20 pcs.
264 203 000	Perfit Neckless: size 3 / cartons of 20 pcs.
264 204 000	Perfit Short: size 4 / cartons of 20 pcs.
264 205 000	Perfit Regular: size 5 / cartons of 20 pcs.
264 206 000	Perfit Tall: size 6 / cartons of 20 pcs.
264 208 000	Carry bag for collars, fits 6 pieces / cartons of 20 pcs.