

Ambu® Res-Cue Pump™



Key Benefits

- Compact and lightweight – small and handy unit – takes up very little space
- Easy to squeeze – requires extremely little force
- Adjustable vacuum (100%, adult or 50%, children) – wide range of application
- Fast and powerful – offering superior performance
- No regular maintenance – low cost operation
- No special storage conditions – stores easily in most emergency kits
- Overfill protection – prevents contamination inside the pump/handle and eliminates need for flush-through procedures and prevents inconvenient exhaust splashing of liquids
- Disposable collection container – no cleaning/sterilisation required

Ambu® Res-Cue Pump™

Power in the palm of your hand

A cost effective alternative to electrical and foot powered devices, the Ambu ResCue Pump is a portable, hand-operated emergency suction pump.

The Ambu Res-Cue Pump is the most versatile hand powered emergency suction pump. A cost effective alternative, offering convenience and performance and a wide range of features/benefits.

Features

- adjustable vacuum.
- Overfill-Protection System prevents contamination of the main unit and averts any inconvenient splashing of fluids.
- very easy to squeeze, requiring little force to obtain effective suction.
- lightweight and robust pump utilises a unique disposable collection container which eliminates the need for cleaning and maintenance.
- Each disposable collection container is supplied with two catheters for adults and paediatrics, ensuring you are prepared for any eventuality.

Specifications

Product	
Vacuum max.	450 mm Hg
Free airflow:	Peak >20l/min
Disposable collection container volume:	300 ml
Stroke adjustment	100% or 50% vacuum
Dimensions	185 x 64 x 168 mm
Weight	230 g

Ordering specifications

Item no.	Description
276 000 001	Ambu Res-Cue Pump/handle complete Includes: Disposable collection container with verfill protection and cap, Rigid suction tubing and adapter – 6 mm, Rigid suction tubing and adapter – 13 mm, 6 mm adapter – for suction catheter and Discard bag
276 000 010	Disposable collection container set Includes: One disposable collection container with overfill protection and cap 2 rigid suction tubings and adapters one 6 mm, one 13 mm and Discard bag
276 000 011	Carry case Can hold the Res-Cue Pump, 2 disposable collection containers and 4 rigid suction tubings



Disposable collection container, suction tubings and adapters

