

# Ambu® TwinPump

	SPARE PARTS	CAT. NO.
1	Suction tube, 1350 mm with combination suction tip and angle Connector	239 000 509
2	Pump chassis, complete with valve diaphragms incl. O-rings and piston rings. 1000 ml container	240 000 501
3	Pump chassis, complete with valve diaphragms incl. O-rings and piston rings. 600 ml container	239 000 501
4	Bottom cover	239 000 514
5	Pedal with retaining springs	239 000 508
6	Piston ring	239 000 516
7	O-ring	239 000 518
8	Combination suction tip	239 000 510
9	Valve diaphragm, separate	239 000 519
10	Manifold pipe	239 000 517
11	Angle connector	239 000 515
12	Cylinder with drawlink and valve diaphragm	239 000 504

