

DATASHEET

Features replaceable skins and latex veins that actually roll when palpated. A realistic flashback confirms proper needle placement, and the shoulder can be used for intramuscular injection training.

For use with all GERi™ and KERi™ manikins.

GERi KERi Manikin IV Arm



KEY BENEFITS

- I.V Arm
- Intramuscular injection training
- Replaceable skins and latex veins
- For use with all GERi™ and KERi™ manikins

ORDERING SPECIFICATIONS

Item no.

LF04080

IV Arm

Replaceable skins and latex veins that actually roll when palpated

Ambu

📍 Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com