

# Ambu® aScope™ 4 RhinoLaryngo Slim

The single-use flexible rhinolaryngoscope



## Key Benefits

- Always ready avoiding any delays
- Improves workflow with no hassle of cleaning and repairs
- High-quality imaging on compact, portable Ambu® aView™ monitor
- Small outer diameter helps minimize patient discomfort during procedures
- Sterile, adding an extra measure of patient safety

## The hassle-free solution that improves your workflow - day and night

From the outpatient clinic to inpatient consult procedures, the single-use, flexible aScope 4 RhinoLaryngo video endoscope can improve your workflow.

### **An excellent choice for nasal endoscopy and laryngoscopy**

The aScope 4 RhinoLaryngo Slim endoscope is designed with a small 3.0mm outer diameter to help minimize patient discomfort during rhinolaryngoscopy procedures. The precise tip motion and high-bending angles allow you to easily manoeuvre the rhinolaryngoscope in the upper airway. The solution includes our compact and high-resolution monitor, Ambu® aView™. When connected to aView, the aScope 4 RhinoLaryngo Slim delivers clear, sharp images allowing you to confidently identify anatomical structures.

And being single-use and sterile, it requires no cleaning, eliminating any delays waiting for a clean scope.

### **Day, night, on weekdays or the weekend**


The aScope 4 RhinoLaryngo endoscope is a flexible, portable, quality solution that can help you be more productive:

- No cost or hassle associated with cleaning and repairs
- Simply take it with you, use it and then discard
- High-quality video imaging that doesn't compromise on the portability you require for consultations outside the clinic
- A brand-new endoscope with the same high quality at every procedure

You can always count on cost transparency with the aScope 4 RhinoLaryngo endoscope.

# Specifications

Optical system	
Field of View	85°
Direction of View	0° (forward viewing)
Depth of Field	6-50 mm
Illumination method	LED
Insertion portion	
Bending section	130° up, 130° down
Insertion cord diameter	3.0 mm (0.12")
Distal end diameter	3.5 mm (0.14")
Maximum diameter of insertion portion	3.5 mm (0.14")
Working length	300 mm (11.8")
Operating environment, storage & transportation	
Temperature - transportation	10 ~ 40° C (50 to 104° F)
Temperature - operating environment	10 ~ 40° C (50 to 104° F)
Temperature - storage	10 ~ 25° C (50 to 77° F)
Relative humidity	30 ~ 85%
Atmospheric pressure	80 ~ 109 kPa
Sterilisation	
Method of sterilisation	ETO




Laryngoscopy performed with aScope 4 RhinoLaryngo Slim and aView



Ambu® aScope™ 4 RhinoLaryngo Slim and Ambu® aView™

## Ordering specifications\*

Item no.	Quantity
510001000 - aScope 4 RhinoLaryngo Slim	5



**Ambu A/S**  
 Baltorpbakken 13  
 DK - 2750 Ballerup  
 Denmark  
 Tel. +45 7225 2000  
 Fax +45 7225 2053  
 www.ambu.com

\* aScope 4 RhinoLaryngo Slim is not available for sale in all countries. Please contact your local Ambu sales representative.