


DESIGNED TO SAVE SPACE AND REDUCE TASK COMPLEXITY

The Ambu airway visualization platform

The Ambu® aView™ 2 and aBox™ 2 endoscopy systems act as the control centre for SureSight and other key Ambu airway visualization devices, all of which are designed to work together.**



**VivaSight and aScope 4 RhinoLaryngo are compatible with aView 2 Advance but not with aBox 2.

Ambu® SureSight™ Connect

Item code	Product	Quantity
390510000	Ambu® SureSight™ Connect	1
390330000	Ambu® SureSight™ Blade MAC Size 3	box of 10
390340000	Ambu® SureSight™ Blade MAC Size 4	box of 10
390230000	Ambu® SureSight™ Blade Hyperangulated Size 3	box of 10
390240000	Ambu® SureSight™ Blade Hyperangulated Size 4	box of 10
390530000	Ambu® SureSight™ Blade Hyperangulated Size 3 Channeled	box of 10

Ambu® aScope™ 5 Broncho/Ambu® aScope™ 5 Broncho Sampler Set

Item code	Product	Quantity
621001000	Ambu® aScope™ 5 Broncho HD 5.0/2.2	5
622001000	Ambu® aScope™ 5 Broncho HD 5.6/2.8	5
620001000	Ambu® aScope™ 5 Broncho 4.2/2.2	5
624001000	Ambu® aScope™ 5 Broncho 2.7/1.2	5
621002000	Ambu® aScope™ 5 Broncho HD 5.0/2.2 Sampler Set	5
622002000	Ambu® aScope™ 5 Broncho HD 5.6/2.8 Sampler Set	5
620002000	Ambu® aScope™ 5 Broncho 4.2/2.2 Sampler Set	5

Ambu® aScope™ 4 Broncho/Ambu® aScope™ 4 Broncho Sampler Set

Item code	Product	Quantity
476 001 000	Ambu® aScope™ 4 Broncho Slim	5
477 001 000	Ambu® aScope™ 4 Broncho Regular	5
478 001 000	Ambu® aScope™ 4 Broncho Large	5
477 501 000	Ambu® aScope™ 4 Broncho Regular Sampler Set	5
478 501 000	Ambu® aScope™ 4 Broncho Large Sampler Set	5

Ambu® VivaSight™ 2

Item code	Product	Quantity
412351000	Ambu® VivaSight™ 2 DLT Kit 35 Fr - LEFT	1
412371000	Ambu® VivaSight™ 2 DLT Kit 37 Fr - LEFT	1
412391000	Ambu® VivaSight™ 2 DLT Kit 39 Fr - LEFT	1
412411000	Ambu® VivaSight™ 2 DLT Kit 41 Fr - LEFT	1
412701000	Ambu® VivaSight™ 2 SLT Kit - 7.0 mm	1
412751000	Ambu® VivaSight™ 2 SLT Kit - 7.5 mm	1
412801000	Ambu® VivaSight™ 2 SLT Kit - 8.0 mm	1
412709000	Ambu® VivaSight™ 2 SLT + EBB Kit - 7.0 mm, 9 Fr	1
412759000	Ambu® VivaSight™ 2 SLT + EBB Kit - 7.5 mm, 9 Fr	1
412809000	Ambu® VivaSight™ 2 SLT + EBB Kit - 8.0 mm, 9 Fr	1
412900000	Ambu® VivaSight™ Endobronchial Blocker Tube	1

Ambu® endoscopy systems

Item code	Product	Quantity
405 011 000	Ambu® aView™ 2 Advance Displaying and Processing unit	1
505001000	Ambu® aBox™ 2 Displaying and Processing Unit	1
HA.3000.991	Ambu® aCart™ Compact Vario Mount	1

Ambu® aScope™ 4 RhinoLaryngo

Item code	Product	Quantity
510 001 000	Ambu® aScope™ 4 RhinoLaryngo Slim	5
512 001 000	Ambu® aScope™ 4 RhinoLaryngo Intervention	5



WANT TO TRY IT?

Your local Ambu rep will call you to arrange a convenient time to show you the product or answer any questions.

BOOK A MEETING NOW



Ambu A/S
Baltorpbakken 13
2750 Ballerup Denmark
Tel.: +45 7225 2000
ambu.com

PUB-002042



Ambu



Ambu® SureSight™ Connect

MAXIMISE YOUR FOCUS ON THE PATIENT

During intubation and airway visualization

Studies in the US have shown that a cluttered workspace and a lack of timely access to familiar tools can increase the time anaesthetists spend switching tasks. This in turn takes the focus away from the patient and increases the risk of errors.¹

Ambu® SureSight™ Connect is a video laryngoscope, which consolidates our trusted airway visualization platform.

It provides a comprehensive range of VL blades for both routine and difficult intubations, and it can be used for emergent or planned procedures.

HELPS YOU MAXIMISE YOUR FOCUS ON THE PATIENT

SureSight gives you access to familiar tools that are designed to simplify time consuming tasks, so you can continue to establish the airway with confidence.

- Use two of our airway visualization devices on the same screen simultaneously
- Familiar, intuitive and integrated tools reduce task complexity
- Connectivity streamlines the sharing and transferring of patient imaging data

¹Jurewicz, et al. Observational study of anaesthesia workflow to evaluate physical workspace design and layout, 2021
Gurses, et al. Identifying and categorising patient safety hazards in cardiovascular operating rooms using an interdisciplinary approach: a multisite study, 2012



EVERYTHING WITHIN REACH AND READY TO GO

Familiar, intuitive, integrated

The Ambu® SureSight™ Connect video laryngoscope is part of the only airway visualization platform that features a one-lung ventilation solution and includes all of the devices recommended for intubation and airway assessment.*

In addition to the video laryngoscope, our trusted single-use aScope bronchoscopes, rhinolaryngoscopes, and VivaSight solutions all offer plug and play compatibility with the Ambu endoscopy systems.**

*Airway visualization equipment recommended to be available for intubation and evaluation: video laryngoscopes, bronchoscopes, rhinolaryngoscopes (UK Difficult Airway Society; American Society for Anesthesiologists; Australian and New Zealand College of Anaesthetists & Faculty of Pain Medicine). The Ambu platform also includes VivaSight DLT, which offers continuous visualization for OLV procedures.

**VivaSight and our rhinolaryngoscopes are compatible with aView 2 Advance but not with aBox 2.



Continuous visualization from two separate devices can reduce task complexity and support your work with airway management techniques and difficult airway algorithms.

USE TWO DEVICES ON THE SAME SCREEN SIMULTANEOUSLY

With the dual view feature in aView 2 Advance or aBox 2**, you are able to simultaneously use two of our airway visualization devices, without needing separate displaying units. For example, you could use the SureSight video laryngoscope with aScope 4 Broncho or another airway visualization device.

**VivaSight and our rhinolaryngoscopes are compatible with aView 2 Advance but not with aBox 2.



SIMPLE TO SHARE AND TRANSFER PATIENT DATA

With the compatible Ambu endoscopy system, you can easily view patient data before the procedure and transfer your newly recorded patient imaging to the PACS after the procedure. If PACS is not available in your department, you can transfer all patient imaging and procedure information files to your hospital's local network drive (SMB).

The connectivity of the Ambu® aView™ 2 Advance and aBox™ 2 endoscopy systems streamlines the task of sharing and transferring patient imaging data. No additional personnel or equipment are needed.