



Ambu

Ambu® SureSight™ Connect

MAXIMISE YOUR FOCUS ON THE PATIENT

During intubation and airway visualization

ONE SOLUTION FOR ALL YOUR AIRWAY VISUALIZATION NEEDS

The Ambu airway visualization platform is designed to save space and reduce task complexity. It is an integrated platform that puts your key airway visualization devices tools within reach and ready to use.

The Ambu airway visualization platform

The only airway visualization platform that features a one-lung ventilation solution and includes all of the devices recommended for intubation and airway assessment.*

*Recommended by the UK Difficult Airway Society; American Society for Anesthesiologists; Australian and New Zealand College of Anaesthetists & Faculty of Pain Medicine



Together with the visualization devices, Ambu laryngeal masks, face masks and resuscitators comprise a comprehensive airway management solution.



Ambu® aView™ 2 Advance Ambu® aBox™ 2

Endoscopy systems - Your control centre for key devices

Ambu® SureSight™ Connect video laryngoscope

Intubation - Maximise your focus on the patient



Ambu® aScope™ 5 Broncho

Bronchoscopy - High performance and flexibility



Ambu® aScope™ 5 Broncho Sampler Sets

A unique integrated BAL/BW solution



Ambu® aScope™ 4 Broncho

Bronchoscopy - Portability and essential performance

Ambu® aScope™ 4 Broncho Sampler Sets

A unique integrated BAL/BW solution

Ambu® aScope™ 4 RhinoLaryngo

Laryngoscopy - For full examination of the upper airway



Ambu® VivaSight™ 2 OLV Solution

OLV procedures - Continuous visualization during your OLV procedure



The Ambu® SureSight™ Connect video laryngoscope solution

Ambu® aCart™ Compact

Everything within reach and ready to go. A compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices



Comprehensive blade portfolio

VL blades for both routine and difficult intubations as well as emergent or planned procedures.

- Mac 3
- Mac 4
- Hyper 3
- Hyper 3C
- Hyper 4

Ambu® aView™ 2 Advance Ambu® aBox™ 2 (not shown here)

Plug and play system that serves as the control centre for key Ambu airway visualization devices

Ambu® SureSight™ Connect

An intuitive, ergonomic and reliable video laryngoscope

- Audible click ensures that the blade is properly attached
- Ergonomic design for secure grip and optimal control
- Anti-fog lens for a clear and consistent image
- Easy-to-reach video/camera button

MAXIMISE YOUR FOCUS ON THE PATIENT

The Ambu® SureSight™ Connect video laryngoscope is part of our trusted airway visualization platform. Having one uniform, intuitive and familiar system helps you stay focused on the patient as you switch between devices and move through the airway management algorithm.

One system, multiple visualization devices

With the Dual View feature in the Ambu endoscopy system, you are able to simultaneously use two of our airway visualization devices without needing separate screens for difficult airways.

Streamlines work with patient imaging data

The Ambu endoscopy system offers connectivity with Electronic Health Records (EHR) to simplify the important task of sharing and transferring patient data.

[Learn more about the Ambu intubation solution](#)

The Ambu® aScope™ Broncho solution

The Ambu® aScope™ Broncho portfolio provides two single-use endoscopy solutions with a variety of sizes available in storage whenever you need them.

Ambu® aView™ 2 Advance Ambu® aBox™ 2 (not shown here)

Plug and play system that serves as the control centre for key Ambu airway visualization devices

Ambu® aCart™ Compact

Everything within reach and ready to go. A compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices



Ambu® aScope™ 5 Broncho

- 2.7/1.2 An ultrathin bronchoscope suitable for ET tubes as small as 3.5 mm
- 4.2/2.2 A regular-sized working channel in a thin bronchoscope
- 5.0/2.2 and 5.6/2.8 Diagnostic and therapeutic scopes for advanced bronchoscopies*



Ambu® aScope™ 4 Broncho

- 3.8/1.2 From difficult intubation to placement and control of DLTs and bronchial blockers in the OR
- 5.0/2.2 From intubation and training to bedside bronchoscopy, including PDT procedures
- 5.8/2.2 With a suctioning channel of 2.8 mm, it can be used for secretion management



Ambu® aScope™ Broncho Sampler Sets - A closed-loop single-use system

- Ambu® aScope™ 4 Broncho Sampler Set
- Ambu® aScope™ 5 Broncho Sampler Set

THE WORLD'S MOST TRUSTED SINGLE-USE BRONCHOSCOPY SOLUTION

It has been tried and tested in millions of cases at thousands of hospitals around the world. It also provides you with a unique, integrated BAL/BW sampling set.

aScope 4 Broncho: Portable, lightweight, easy to set up

Offers the portability and essential performance needed for a variety of procedures in the OR to routine bronchoscopy in the ICU

aScope 5 Broncho: A new standard for single-use bronchoscopes

Combines the high level of manoeuvrability and imaging needed or advanced diagnostic and therapeutic procedures with the sterility and efficiency of the single-use concept

aScope Sampler Set: Designed for superior lavage

Reduces the risk of sample contamination and loss common with Lukens traps, protecting the sample from start to finish

[Learn more about the Ambu bronchoscopy solutions](#)

The Ambu® VivaSight™ 2 OLV solution

Ambu® aCart™ Compact

Everything within reach and ready to go. A compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices

Ambu® VivaSight™ 2 DLT

Built-in camera at the tip of the tracheal lumen enables continuous visualization

Ambu® VivaSight™ 2 SLT*

Integrated camera enables real-time monitoring of correct tube placement and continuous visualization of endobronchial blocker position

VivaSight Endobronchial Blocker Tube*

Designed for One Lung Ventilation when used in conjunction with VivaSight 2 SLT

Ambu® aView™ 2 Advance

Plug and play system that serves as the control centre for key Ambu airway visualization devices

Ambu® aScope™ 4 and 5 Broncho

A variety of different sizes are always available to handle the various situations that can arise during OLV procedures



CONTINUOUS VISUALIZATION INCREASES PATIENT SAFETY AND IMPROVES WORKFLOW

The Ambu® VivaSight™ 2 portfolio offers a safe and cost-effective way to perform One Lung Ventilation (OLV). Although similar to conventional double lumen tube (DLT) and single-lumen tube (SLT) solutions, VivaSight 2 has a built-in camera at the tip of the tracheal lumen.

This enables you to immediately detect and correct tube malpositioning and dislocation without compromising ventilation of the patient.

No need to set up a second system

VivaSight 2 significantly reduces the need to use a bronchoscope. In cases where you do need a bronchoscope, however, you can take advantage of Dual View with Ambu® aScope™ 4 or 5 Broncho and aView 2 Advance (Gen 2). There's no need to set up second system.

[Learn more about the Ambu VivaSight 2 OLV solution](#)

*Not available in the US

The Ambu® aScope™ RhinoLaryngo solution


The Ambu® aScope™ RhinoLaryngo portfolio provides two single-use endoscopy solutions available in storage whenever you need them.

Ambu® aView™ 2 Advance

Plug and play system that serves as the control centre for key Ambu airway visualization devices

Ambu® aCart™ Compact

Everything within reach and ready to go. A compact, mobile workstation for the SureSight Connect video laryngoscope solution and your other key devices



Ambu® aScope™ 4 RhinoLaryngo

- 3.0 outer diameter: An excellent choice for nasal endoscopies, laryngoscopies and swallow assessments
- 5.0/2.0 An excellent choice for therapeutic procedures

HELPS YOU STREAMLINE YOUR WORKFLOW AND PERFORM PROCEDURES CONFIDENTLY

The Ambu® aScope™ 4 RhinoLaryngo solution is a single-use endoscopy solution that offers excellent imaging, is always available and easy to take with you.

Excellent imaging with integrated video and audio

- Easily identify anatomical structures with high endoscope image quality and the Ambu full-HD endoscopy system
- Dictate notes as you perform the procedure with synchronised audio and video
- Review recordings of fast movements more carefully with slow-motion video playback
- Isolate and save key frames for reporting, training and to get second opinions

[Learn more about the Ambu RhinoLaryngoscope solution](#)

THE AMBU AIRWAY VISUALIZATION PLATFORM

One single-use platform; multiple long-term benefits

The Ambu endoscopy systems are used across many clinical areas. This means you can use the familiar Ambu endoscopy system already present in the department when you are called to establish an airway.

Having one high-quality platform designed to work together means:

- A familiar and intuitive system
- One trustworthy supply partner to simplify procurement and provide ongoing support
- Access to a portfolio of single-use endoscopes and OLV devices, which can reduce costs associated with capital investments, reprocessing and repairs

ORDERING NUMBERS

Ambu® SureSight™ Connect

| Item code | Product | Quantity |
|-----------|--|-----------|
| 390510000 | Ambu® SureSight™ Connect | 1 |
| 390330000 | Ambu® SureSight™ Blade MAC Size 3 | box of 10 |
| 390340000 | Ambu® SureSight™ Blade MAC Size 4 | box of 10 |
| 390230000 | Ambu® SureSight™ Blade Hyperangulated Size 3 | box of 10 |
| 390240000 | Ambu® SureSight™ Blade Hyperangulated Size 4 | box of 10 |
| 390530000 | Ambu® SureSight™ Blade Hyperangulated Size 3 Channeled | box of 10 |

Ambu® aScope™ 5 Broncho/Ambu® aScope™ 5 Broncho Sampler Set

| Item code | Product | Quantity |
|-----------|--|----------|
| 621001000 | Ambu® aScope™ 5 Broncho HD 5.0/2.2 | 5 |
| 622001000 | Ambu® aScope™ 5 Broncho HD 5.6/2.8 | 5 |
| 620001000 | Ambu® aScope™ 5 Broncho 4.2/2.2 | 5 |
| 624001000 | Ambu® aScope™ 5 Broncho 2.7/1.2 | 5 |
| 621002000 | Ambu® aScope™ 5 Broncho HD 5.0/2.2 Sampler Set | 5 |
| 622002000 | Ambu® aScope™ 5 Broncho HD 5.6/2.8 Sampler Set | 5 |
| 620002000 | Ambu® aScope™ 5 Broncho 4.2/2.2 Sampler Set | 5 |

Ambu® aScope™ 4 Broncho/Ambu® aScope™ 4 Broncho Sampler Set

| Item code | Product | Quantity |
|-------------|---|----------|
| 476 001 000 | Ambu® aScope™ 4 Broncho Slim | 5 |
| 477 001 000 | Ambu® aScope™ 4 Broncho Regular | 5 |
| 478 001 000 | Ambu® aScope™ 4 Broncho Large | 5 |
| 477 501 000 | Ambu® aScope™ 4 Broncho Regular Sampler Set | 5 |
| 478 501 000 | Ambu® aScope™ 4 Broncho Large Sampler Set | 5 |

Ambu® VivaSight™ 2

| Item code | Product | Quantity |
|-----------|---|----------|
| 412351000 | Ambu® VivaSight™ 2 DLT Kit 35 Fr - LEFT | 1 |
| 412371000 | Ambu® VivaSight™ 2 DLT Kit 37 Fr - LEFT | 1 |
| 412391000 | Ambu® VivaSight™ 2 DLT Kit 39 Fr - LEFT | 1 |
| 412411000 | Ambu® VivaSight™ 2 DLT Kit 41 Fr - LEFT | 1 |
| 412701000 | Ambu® VivaSight™ 2 SLT Kit - 7.0 mm | 1 |
| 412751000 | Ambu® VivaSight™ 2 SLT Kit - 7.5 mm | 1 |
| 412801000 | Ambu® VivaSight™ 2 SLT Kit - 8.0 mm | 1 |
| 412709000 | Ambu® VivaSight™ 2 SLT + EBB Kit - 7.0 mm, 9 Fr | 1 |
| 412759000 | Ambu® VivaSight™ 2 SLT + EBB Kit - 7.5 mm, 9 Fr | 1 |
| 412809000 | Ambu® VivaSight™ 2 SLT + EBB Kit - 8.0 mm, 9 Fr | 1 |
| 412900000 | Ambu® VivaSight™ Endobronchial Blocker Tube | 1 |

Ambu® endoscopy systems

| Item code | Product | Quantity |
|-------------|---|----------|
| 405 011 000 | Ambu® aView™ 2 Advance Displaying and Processing unit | 1 |
| 505001000 | Ambu® aBox™ 2 Displaying and Processing Unit | 1 |
| HA.3000.991 | Ambu® aCart™ Compact Vario Mount | 1 |

Ambu® aScope™ 4 RhinoLaryngo

| Item code | Product | Quantity |
|-------------|---|----------|
| 510 001 000 | Ambu® aScope™ 4 RhinoLaryngo Slim | 5 |
| 512 001 000 | Ambu® aScope™ 4 RhinoLaryngo Intervention | 5 |

THE AMBU COMMITMENT

Ambu is on a mission to deliver airway management solutions that help you care for patients and save lives. At the same time, as a single-use medical device company, we are acutely aware of our obligation to work towards more sustainable solutions.

Taking leaps forward in sustainability
Advancing in the area of sustainability is one of the four strategic pillars of our company. This strategy focuses on circular design principles that help us lower our products' carbon footprint and activities throughout our value chain that help us work towards our net zero emissions goal.



One concrete example of this is our pioneering use of bioplastic material, which is a mix of fossil-based and 2nd gen. bio feedstock, in our endoscope handles, laryngeal mask cuff protector packaging and other products.

Our aim is that by the end of 2024, all our currently marketed endoscopes will be made with bioplastic material in the handles. In the future, we plan to introduce the material in other parts of our endoscopes, too. In addition, SureSight is phthalate and PVC-free.

Contact us for more information

Learn more about the Ambu airway visualization platform, including solutions for intubation, bronchoscopy and One Lung Ventilation.

[Learn more about the Ambu Airway Visualization Platform](#)