







ENDOSCOPY SYSTEMS

ENDOSCOPE NAME	CONNECTOR CONFIGURATION	Ambu® aView™ 2 Advance 1 st Gen	Ambu® aView™ 2 Advance 2 nd Gen	Ambu® aView™ 2 Advance 2 nd Gen	Ambu® aBox™ 2
					
VivaSight™ 2 SLT/DLT		✓	✓	✓	
aScope™ 4 RhinoLaryngo family		✓	✓	✓	
aScope™ 4 Broncho family		✓	✓	✓	
aScope™ 4 Cysto		✓	✓	✓	
aScope™ 5 Cysto HD			✓	✓	✓
aScope™ 5 Uretero			✓		✓
aScope™ 5 Broncho HD 5.0/2.2 & 5.6/2.8			✓	✓	✓
aScope™ 5 Broncho 2.7/1.2 & 4.2/2.2			✓		✓
SureSight™ Connect			✓	✓	✓
aScope™ Gastro/Gastro Large					✓
aScope™ Colon					✓
aScope™ Duodeno 2					✓