

Ambu® ResCue Mask



Key Benefits

- Transparent dome allows rescuer to visually check for patient's lip colour and vomitus
- Oxygen inlet (optional) facilitates increased oxygen delivery to patient
- Made of durable plastics for easy cleaning and long product life
- Private labelling available
- Available in hardcase, softcase or plastic bag
- Phthalate-free material

Ambu® ResCue Mask

The Ambu ResCue Mask is designed to provide mouth-to-mask ventilation.

As a physical barrier between the rescuer and the patient Ambu ResCue mask is used to give protected rescue breathing for any patient (any age group) and may also be used with a resuscitator.

The Ambu ResCue Mask helps protect rescuers in an emergency situation. It prevents direct contact with patient's mouth, nose and face and helps to overcome hesitation to start resuscitation.

Specifications

Product	
Size	
Dimensions	13.8 x 11 x 4.8 cm
Weight hardcase	112 g
Weight softcase	67 g
Materials	
Mask	Polyvinylchloride
Softcase	Nylon
Hardcase	Polypropylene



Ambu Res-Cue Mask hardcase and softcase.

Ordering specifications

Item no.	Ambu Rescue Mask
000 252 225	Ambu Res-Cue Mask with O ₂ inlet, in red hard case
000 252 226	Ambu Res-Cue Mask without O ₂ inlet, in red hard case
000 252 223	Ambu Res-Cue Mask with O ₂ inlet, in soft case
000 252 224	Ambu Res-Cue Mask without O ₂ inlet, in soft case
000 252 222	Ambu Res-Cue Mask without O ₂ inlet, in plastic bag
000 252 253	One way valve with bacterial filter
000 252 123	Ambu Res-Cue Mask – Adult and Infant mask in soft case without O ₂ Inlet
Item no.	Ambu Rescue Mask Private Label
000 252 235	Ambu Res-Cue Mask private label hard case with O ₂ inlet
000 252 236	Ambu Res-Cue Mask private label hard case without O ₂ inlet

*ResCue Mask with Oxygen Inlet is Rx Only in the USA



Ambu Res-Cue Mask with Adult and Infant mask in a softcase.

Logo
Organisation, company or sponsor logo. Maximum logo size: 50 x 40 mm
Hardcase
Available in all pantone colours
Print
Four colour prints are available on either a white or transparent sticker
1 colour print if logo is printed on the case itself
Order
Minimum 250 pieces



Ambu Res-Cue Mask private label.