

# Ambu® ResCue Mask



## Key Benefits

- The ResCue Mask is transparent to allow visual check of patient's condition
- Head strap for fixation and secure positioning. For use on both adults and children
- Reusable facemask that is easy to clean, which extends the product life
- Private labelling available
- Available in softcase

## Ambu® ResCue Mask

The Ambu ResCue Mask is designed to provide protected rescue breathing for any patient (any age group). In CPR, mouth-to-mouth has been the traditional way of ventilation in an emergency situation. With the Ambu ResCue Mask first responders can provide immediate respiratory support without having any concerns about infectious diseases and the inherent unhygienic nature of the practice thus removing any natural reluctance to carry it out.

### Features

The Ambu ResCue Mask helps protect rescuers in an emergency situation. It prevents direct contact with patient's mouth, nose and face and helps overcome hesitation to start resuscitation.

The mask is available with and without oxygen inlet. Thus an oxygen tube can be attached to the oxygen inlet for supplemental oxygen. The one-way filter of the mask can be removed in order to use the mask with a manual resuscitator.

The Ambu ResCue Mask comes in an original modern shell-shaped softcase.

# Specifications

<b>Product</b>	
<b>Size</b>	
<b>Dimensions</b>	138 mm x 113 mm x 50 mm
<b>Weight</b>	Approximately 110g
<b>Materials</b>	
<b>Mask</b>	Polyvinylchloride
<b>Softcase</b>	Polyurethane, Ethylene Vinyl Acetate and Micro Fabric



The mask is packed inside the softcase with gloves and alcohol pad.

## Ordering specifications

<b>Item no.</b>	<b>Ambu ResCue Mask</b>
252 301 004	Ambu ResCue Mask – Adult mask with oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase
252 301 104	Ambu ResCue Mask – Adult mask without oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase
<b>Item no.</b>	<b>Ambu ResCue Mask Private Label</b>
252 320 004	Ambu ResCue Mask – Adult mask with oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase
252 320 104	Ambu ResCue Mask – Adult mask without oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase

\*ResCue Mask with Oxygen Inlet is Rx Only in the USA

<b>Logo</b>
Organisation, company or sponsor logo. Maximum logo size: 70 x 40 mm
<b>Order</b>
Minimum 12000 pieces