

# Ambu® Airway Man W + I



## Key Benefits

- Wireless connections
- Ambu intubation head
- The rigidity of the chest is adjustable
- Airways open only when head is correctly extended
- Jaw thrust possible
- PC compatible
- Training mat integrated in the carrying bag

## Ambu® Airway Man W + I

The Ambu Airway Man is a highly developed instruction- and training manikin that provides a remarkably lifelike representation of the human anatomy, particularly those features important to training in modern resuscitation.

### Connection free

The Airway Man does not need cables or wires, which means it can be placed away from computers and power supplies. This gives a great amount of freedom while training and because there are no wires or cables lying around it gives an even more realistic scenario.

### Mechanical Instrument

#### Version W:

The Airway Man W (Wireless) models are equipped with a mechanical monitoring instrument as well as digital sensors showing effectiveness of the resuscitation.

It also shows results achieved for insufflations volume and depth of external chest compression, as well as it indicates any stomach inflation and incorrect hand positioning.

Airway Man W models can be connected to a computer through USB or WLAN using the Ambu CPR-Software for a more detailed analysis of the performed resuscitation.

#### Version I:

The Airway Man I (Instrument) models are equipped with a mechanical monitoring instrument showing effectiveness of the resuscitation, including results achieved for ventilation volume and depth of external chest compression. Moreover it indicates any stomach inflation and incorrect hand positioning.

## Specifications

Product	
Torso with carrying case	13 kg
Total length: Torso	80 cm

## Materials

Product	Material
Basic unit	Polyethylene
Instrument part	ABS plastic
Skull	PVC, hard
Chest, face	PVC, soft



Ambu Airway Man with an AuraOnce

## Ordering specifications

Item no.	Description
A284 224 000	Ambu Airway Man I, Torso
A284 327 000	Ambu Airway Man W, Torso with CPR-Software
234 000 706	Upgrade Kit - Full Body
234 000 704	Upgrade Kit - Arms
Accessories	
234 000 790	Jogging suit jacket
234 000 792	Jogging suit trousers
234 000 786	Legs with trousers
234 000 794	Carrying bag



Bag for Ambu Airway Man