



Ambu® aScope™ 4 RhinoLaryngo

**The hassle-free solution
that improves your workflow
-day and night**

A flexible, portable, quality solution

From the outpatient clinic to inpatient consult procedures, the single-use, flexible aScope 4 RhinoLaryngo video endoscope can improve your workflow. It can help you be more productive—day, night, on weekdays or the weekend.

Just grab it and go

Always ready when you are, the single-use aScope 4 RhinoLaryngo endoscope requires no reprocessing. Simply use it and discard. No more delays waiting for a clean scope; no cost and hassle associated with cleaning and repairs.

The solution includes our compact, high-resolution monitor, Ambu® aView™, which is intuitive to use and fits easily into small spaces.

It is lightweight and simple to setup, so you can take it with you wherever you need to go in the hospital or clinic.



High-quality imaging you can rely on

The aScope 4 RhinoLaryngo solution brings you high-quality video imaging on an external monitor without compromising on the portability you require for consultations outside the clinic.

When connected to the aView monitor, the aScope 4 RhinoLaryngo gives you clear visibility of the upper airway, and the sharp images and adaptive light allow you to confidently identify anatomical structures.

Navigation in the upper airway is made easy thanks to its high-bending angles and precise tip motion. And since it is single-use, you get a brand-new endoscope with the same high quality at every procedure.





Brings you closer to your patients and peers

The aView monitor makes it easier to focus on and involve your patients in the procedure. Compact and lightweight, you can simply manoeuvre it into a comfortable position close to you and the patient. In addition, having a separate monitor promotes a better working posture compared to traditional fibrescopes with eyepieces. All in all, the complete solution makes for a more pleasant examination

experience for you and your patients. It is also simpler to train others and share knowledge. Your peers can follow the procedure on the monitor, and you can record and store images and videos of the procedures for later viewing and documentation.

Lightweight, ergonomic design

No matter how many procedures you do each day, the aScope 4 RhinoLaryngo endoscope will feel comfortable in your hands. The lightweight design, with its ergonomic, rippled-surface handle, makes it easy for you to manoeuvre during examinations.


Guaranteed sterility, single-use efficiency

When it comes to disinfection techniques for reusable endoscopes, guidelines within ENT are becoming increasingly strict. The single-use, sterile aScope 4 RhinoLaryngo endoscope allows you to adhere to these guidelines and improves your workflow. Regardless of the degree of sterility needed for your upper airway procedures, a single-use, sterile scope simply adds an extra measure of patient safety.

The simple way to improve productivity and safety

Being sterile, the aScope 4 RhinoLaryngo endoscope eliminates any risk of device-related cross-infection. You simply use the scope and discard it. There is no need to establish a

new, more complicated workflow or to worry about the costs associated with cleaning and maintaining reusable scopes. You can always count on cost transparency with the aScope 4 RhinoLaryngo endoscope.



A complete solution that improves your everyday workflow

The aScope 4 RhinoLaryngo solution includes the compact, high-resolution aView monitor and a portfolio of two single-use flexible video rhinolaryngoscopes:

● The aScope 4 RhinoLaryngo Intervention*

The aScope 4 RhinoLaryngo Intervention endoscope is designed for therapeutic procedures. It is ready when you are, requires no post-procedural handling and comes at a low investment cost, which makes it particularly suitable for outpatient-clinic biopsies. The Intervention delivers high-quality imaging and easy manoeuvrability.

- **Outer diameter:** 5.0 mm
- **Working channel:** 2.2 mm
- **Length:** 350 mm
- **Bending angles:** 130°/130°

● The aScope 4 RhinoLaryngo Slim*

The aScope 4 RhinoLaryngo Slim endoscope is an excellent choice for nasal endoscopy and laryngoscopy. It delivers high-quality imaging on the aView monitor, and the small outer diameter helps minimize patient discomfort during procedures. The precise tip motion and high-bending angles allow you to easily manoeuvre the endoscope in the upper airway.

- **Outer diameter:** 3.0 mm
- **Working channel:** N/A
- **Length:** 300 mm
- **Bending angles:** 130°/130°



Ambu A/S

Baltorpbakken 13
DK-2750 Ballerup
Denmark

For more information

Visit ambu.com

Trusted innovation from Ambu

Since 1937 Ambu has been bringing breakthrough healthcare solutions to life within anaesthesia, patient monitoring and diagnostics, and emergency care. Today, we are dedicated to improving safety for millions of patients and healthcare professionals around the world through the advancement of single-use devices, including aScope – the world's first single-use flexible video scope.

Our commitment to bringing new ideas, exceptional product performance and superior service to our customers has made Ambu one of the most recognised medical companies in the world. Headquartered near Copenhagen, Denmark, Ambu employs around 2,500 people across Europe, North America and Asia Pacific. For more information, please visit www.ambu.com

*** aScope 4 RhinoLaryngo Intervention and aScope 4 RhinoLaryngo Slim are not available for sale in all countries. Please contact your local Ambu sales representative.**